innoliving

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UPPER ARM AUTOMATIC DIGITAL BLOOD PRESSURE MONITOR

GB

Thank you for purchasing the model INN-006 Upper Arm Automatic Digital Blood Monitor

For a proper use of the product, please read the following instructions carefully and save this manual for possible future consultation.

This instruction manual is intended to assist the user for safe and efficient operation of the automatic digital blood pressure monitor INN-006. The device must be used in accordance with the procedures described in the manual. It is important to read and understand the entire manual, especially the section < **IMPORTANT SAFETY INSTRUCTIONS** >.

This device is intended for the non-invasive measurement of systolic and diastolic arterial blood pressure and pulse rate in adults (age 15 and above).

CAUTION

1. Do not use this device on infants or persons who cannot express their intentions.
2. The device is not suitable for measuring the blood pressure of children. Ask your doctor before using it on elder children.
3. The patient is an intended operator. But persons who suffer from arrhythmia, diabetes, cardiovascular problems or who have had a stroke should consult your doctor before using the device.

PRINCIPLE OPERATION

This device adopts the oscillometric technology with Fuzzy Algorithm to measure the arterial blood pressure and pulse rate. The cuff is wrapped around the arm and automatically inflated by the air pump. The sensor of the device catches weak fluctuation of the pressure in the cuff produced by extension and contraction of the artery of the arm in response to each heartbeat. The amplitude of the pressure waves is measured, converted in millimeters of the mercury column, and is displayed by digital value.

ATTENTION: This device can not provide reasonable accuracy if used or stored in the temperature, humidity or altitude beyond the range stated in the section <SPECIFICATIONS> of this manual.

NEW TECHNOLOGIES USED

Fuzzy Algorithm is the processing algorithm, taking into account the specialty of individual heartbeats, which provides higher accuracy of measurement. Software version: V1.1

IMPORTANT SAFETY INSTRUCTIONS

It is necessary to know that arterial blood pressure is subjected to sharp fluctuations. The level of the arterial blood pressure depends on many factors. Generally arterial blood pressure is lower in summer and higher in winter. Arterial blood pressure changes with atmosphere pressure and is affected considerably by many factors, e.g. physical loads, emotional excitability, stress, meals, etc. Medicines, drinking, smoking affects greatly the level of individual blood pressure. Blood pressure does vary with age and individuals, and it is recommended to write down the readings from blood pressure records daily, then you can check with your doctor to find out what is a "normal blood pressure measurement" for you.

Please read the instruction manual carefully before using this device especially <IMPORTANT SAFETY INSTRUCTIONS>, it can help you use the device correctly and safely! Please keep the instruction manual safe for future use. For specific information about your own blood pressure, consult your physician.

WARNINGS

- Consult your physician if you suffer from illnesses prior to using the device.
- The device is not suitable for persons who have electrical implants.
- If you had a mastectomy (breast amputation) do not use this blood pressure monitor on the arm on the side of the mastectomy.
- Pregnant women should only measure their own blood pressure in consultation with their doctor, since the readings may be changed with pregnancy.

- Do not service or maintain the cuff and the device while in use with patient.
- Do not use this blood pressure monitor on any arm where intravascular access or therapy (such as an intravenous drip or a blood transfusion), or an arteriovenous shunt (A-V shunt) is present. The temporary interference to blood flow by the blood pressure measurement could result in injury.
- Do not use the device with other medical electrical (ME) equipment simultaneously.
- Do not use the device in the area the HF surgical equipment, MRI, or CT scanner exists, or in the oxygen rich environment.
- Do not use a mobile phone or other devices that emit electromagnetic fields, near the device. This may result in incorrect operation of the device.
- Verify (for example, by observation of the limb concerned) that operation of the AUTOMATED SPHYGMOMANOMETER does not result in prolonged impairment of the circulation of the blood of the patient.
- Never use any accessories or parts from other manufacturers. Using such accessories or parts could cause a hazardous situation for the user or damage to the device.
- Do not modify this device without authorization of the

manufacturer.

- The batteries used in this device may present a fire or chemical burn hazard if mistreated. Do not disassemble, heat or incinerate.
- Keep this device away from fire and heat sources to prevent fire or explosion
- Please keep the unit out of reach of infants, children or pets, since inhalation or swallowing of small parts can be dangerous or even fatal.
- Please pay attention that the continuous CUFF pressure due to connection tubing kinking will cause a harmful injury.
- Do not use an extension cord with this device.
- The air tube or the AC adapter cable may cause accidental strangulation in infants.
- Do not put the air tube around your neck - danger of suffocation!
- A device should never be left unattended when plugged in.
- Do not reach for a corded device that has fallen into water. Unplug immediately.
- It is quite normal that two measurements taken in quick succession may produce significantly different results, because too frequent and consecutive measurements could cause disturbances in blood circulation and injuries.
- This medical device is intended for household use only.

CAUTIONS

- Use this device under the right environmental conditions as indicated in this user manual. If not, this could affect the performance, lifetime of the device and measurement results.
- Only use this device for its intended purpose as described in this user manual.
- Do not confuse self-monitoring with self-diagnosis. This device allows you to monitor your blood pressure. Do not begin or end medical treatment based on the measurement results. Always consult your physician for treatment advice.
- Do not take any therapeutic measures on the basis of a self-measurement. Never change prescribed medication without consulting your physician. Consult your physician if you have any questions about your blood pressure.
- In case you will obtain unexpected readings, consult your physician.
- If you are taking medication, consult your physician to determine the most appropriate time to measure your blood pressure.
- Consult the physician if measurement errors occur in children or persons with arrhythmia.
- The pulse display is not suitable for monitoring the frequency of cardiac pacemakers.

- Common arrhythmias (such as atrial or ventricular premature beats or atrial fibrillation) and peripheral artery disease / arteriosclerosis can affect the accuracy of this blood pressure monitor. Please consult your physician how to best use this blood pressure monitor if you suffer from any of these conditions. Blood pressure measurement is not suitable in cases of serious arteriosclerosis (hardening of the arteries).
- The effectiveness of this blood pressure monitor has not been established in pregnant women.
- Always check the device and cuff before you use it. Do not use the device or cuff if one of them is damaged, because this may cause injury.
- This device is not intended for use on extremities other than the arm or for functions other than obtaining a blood pressure measurement.
- Do not attach the cuff on the same arm on which other monitoring medical electrical equipment is attached simultaneously, because this could cause temporary loss of function of those simultaneously-used monitoring medical electrical equipment.
- Never attach the cuff on injured skin, an injured arm or an arm under medical treatment as this can cause further injury.
- Do not forcibly crease the arm cuff or the air tube

excessively.

- Do not press the air tube while taking a measurement.
- Do not use the device in case of existing polyester or nylon material allergies.
- This device is not suitable for continuous monitoring during medical emergencies or operations.
- This device cannot be used with HF (High Frequency) surgical equipment at the same time.
- This device is not washable. Never immerse the device in water and do not rinse it under the tap.
- This device should keep dry to prevent from moisture.
- The equipment is not AP/APG equipment and is not suitable for use in the presence of a flammable anesthetic mixture with air, with oxygen or nitrous.
- To avoid measurement errors, do not use the device near strong electromagnetic fields, radiated interference signal or electrical fast transient/burst signal. For example magnets, radio transmitters, microwave ovens.
- If this device was stored in low temperature, leave it in room temperature for at least 1 hour.
- Repeated measurements with an interval of 5 minutes are recommended, so you can calculate the average to get a more accurate measurement. An interval of 5 minutes can also ensure that the operation of the device does not result in prolonged impairment of the

circulation of the blood.

- Atherosclerosis patients may require longer interval (10-15minutes) as elasticity of patient's vessels decreases significantly with the disease. 10-15minutes interval is also applicable for patients suffering from diabetes for a long period of time.
- Dispose of the device, components and optional accessories according to applicable local regulations. Unlawful disposal may cause environmental pollution.
- Connecting electrical equipment to mso effectively leads to creating a ME system, and can result in a reduced level of safety.

CONTRAINDICATIONS AND WARNINGS FOR USE

As it happens for all blood pressure monitors, some particular conditions might affect the measurement, among the others:

- Heart rhythm disorders
- Kidney disorders
- Very low blood pressure
- Very low blood perfusion
- Patient in shock
- Diabetes
- Abnormalities of the blood vessel
- People with electrical systems like a cardiac pacemaker
- Pregnant women


- Women with preeclamptic pathology
- Patient movement during measurement

These conditions can produce incorrect readings, which represent a danger to your health, as values can be misinterpreted. Always consult your doctor to assess your health appropriately.

CLASSIFICATION

- ME EQUIPMENT not intended for use in an oxygen rich environment or in the presence of flammable mixers.
- Internally powered equipment (without adaptor), Class II equipment (with adaptor).
- Type BF applied part, recognize the cuff as applied part.

PARTS AND COMPONENTS

1. Main Body
2. LCD Display
3. Air Connector
4. Tube Plug
5. Air Hose
6. Cuff
7. D-ring
8. Button M
9. Button 
10. Batteries (Optional)

BATTERY INSTALLATION

1. Open the battery cover and install four 'AA' type batteries into the battery compartment as indicated. Make sure that the polarity is correct;
2. Close the battery compartment cover.
 - Replace the batteries when the replacement indication  appears in the display or nothing after  button is pressed;
 - Batteries in this kit are intended to check work capacity of the device and the life-span of the batteries can be shorter than the recommended time;
 - Use R6, LR6 or AA alkaline batteries, do not use rechargeable batteries;
 - Only same type batteries are allowed to be used together;
 - Replace all batteries simultaneously;
 - If the device is to be unused for a long time, please take out the batteries;
 - Don't leave the worn batteries in the device.

USE THE DEVICE WITH AC POWER ADAPTER

Besides batteries you can use AC power adapter as the power supply. The AC power adapter is optional for this device for sale.

The AC adapter is specified as a part of the blood pressure monitor.

- Insert the AC adapter cord into the jack on the right side of the monitor
- Insert the AC adapter plug into the outlet.
- To remove the AC adapter, disconnect the adapter plug from the AC outlet first and then disconnect the cord from the monitor's jack.

CAUTION

- When using optional AC adapter, the AC adapter must comply with the requirements of standard IEC60601-1.

- To avoid possible damage of the monitor, use only the exclusive AC adapter that can be purchased from authorized dealers. Other adapter may damage the blood pressure monitor.
- The AC adapter is used as an isolating means, the AC adapter plug shall insert into the outlet nearby the operator, make it easy to disconnection the device from the outlet.
- If long time work, remove the plug after the adapter cools, and prevent burns.
- Plug the AC adapter into the appropriate voltage outlet. Do not use in a multi-outlet plug.
- Not to position the blood pressure monitor to make it difficult to operate the disconnection device(adaptor).

NOTE: The monitor is designed not to draw power from the batteries when the AC adapter in use.

CAUTION

- Please keep quiet for 5-10 minutes and avoid eating, drinking, alcohol, smoking, exercising and bathing before taking a measurement. All these factors will influence the measurement result.
- Remove any garment that fits closely to your upper arm.
- Always measure on the same arm (normally left).
- Measurements should be taken regularly at the same time of each day, as the blood pressure varies even during the day.
- Any effort to support the arm during measurement may increase the measured blood pressure.
- Make sure, you are in a comfortable, relaxed position with leg-uncrossed, feet flat on the floor, back and arm supported, middle of the cuff at the level of the right atrium of the heart and do not move or constrict your muscles and talk during measurement. Use a cushion to support your arm if necessary. Keep position in normal use.
- If the arm artery lies lower or higher than the heart, a false reading will be obtained.
- A loose or open cuff causes false readings.
- With repeated measurements, blood accumulates in the arm which can lead to false reading.
- Consecutive blood pressure measurements should be repeated after 1 minute pause or after the arm has been held up in order to allow the accumulated blood to flow away.

Optional AC adapter technical feature

Model: Y55M-0600600

Input: 100-240V 50/60Hz

Output voltage: $6V \pm 5\%$

Output current: 600 mA

Output plug polarity: $\leftarrow \rightarrow$ inner

Insulation degree: II


Fig.1

CORRECT POSTURE FOR MEASUREMENT

1. Sit beside the table and let the table support your arm as you take the measurement. Make sure that the cuff on the upper arm is at approximately the same level as the heart, and that the forearm is extended naturally on the table as in Fig.1.

2. You may lie on your back and take a measurement. Look at the ceiling, keep calm, and don't move your neck or body during the measurement. Make sure that the cuff on the upper arm is at approximately at the same level as the heart as in Fig.2.





Fig.2

ASSEMBLING THE CUFF

1. Insert the edge of the cuff approximately 5 centimeters into the D-ring as shown in figure.


2. Put the cuff on the left upper arm with the tube pointing to the direction of the palm. If measurement on your left arm is difficult, you can use the right arm for measurement. In this case, it is necessary to know that the readings may differ about 5-10mmHg between left arm and right arm.


3. Wrap cuff around your upper arm with the lower

edge of the cuff approximately 2-3 centimeters above the elbow. The mark "ARTERY" must be over the artery of the arm.


4. Press the cuff to make sure that it is attached securely. The cuff should not be too tight or loose is greatly recommended. Two fingers should be easily put in between cuff and upper arm.


5. The mark INDEX on the cuff must point to area NORMAL. This means the cuff size is correct. If mark INDEX points to area beyond area NORMAL, please consult your dealer whether you need another size cuff.

6. Sometimes is difficult to make the cuff regular depending in the shape of the user's upper arm, the cone-shape assembly of cuff is also acceptable.

7. If your clothes restrict the blood circulation of your upper arm, or you roll your sleeve up so as to result in such restriction. Please take of your shirt to get an accurate measurement if necessary.

CARRYNG OUT A MEASUREMENT

1. Assembling the cuff, then insert the tube plug into the air connector. Before the measurement, take 3-5 times deep breathe and relax yourself. Don't talk or move your arm.

2. Press the button  and all symbols will appear on display in 2 seconds as in Fig.1.

Fig. 1

Then two short beeps will sound and “0 mmHg” will appear on the screen. Pump begins to inflate with display showing the reading of pressure. Generally the pressure will reach 190 mmHg, as in Fig. 2.

Fig. 2

3. **NOTE:** The device will inflate to a higher pressure automatically in case the inflation pressure is not sufficient to determine measurement result.

4. The pump stops inflating and the pressure begins to decrease gradually, during which the user's blood pressure and pulse will be calculated as in Fig. 3.

Fig. 3

5. There will be a long beep following the accomplishment of measurement. The air in the cuff will deflate quickly and the blood pressure reading, pulse reading will show in the display as in Fig. 4


Fig. 4

6. If irregular heartbeat was detected during the measurement, LCD display the  icon to remind users of heartbeat irregularity as Fig. 5.

ATTENTION: We recommend contacting your physician if you see the  indicator frequently.


Fig. 5

7. Press the button  to turn off the device. Please rest for at least 3 minutes for another measurement. If the power supply is not switched off and the device keeps unused for 3 minutes, the device will be switched off automatically.

FUNCTION OF MEMORY

Memory recall

1. The device can store 90 sets of readings respectively in M1 and M2, and will automatically calculate the average value of the latest 3 readings for M1 and M2. When the memory is full (90 sets of readings are stored), the oldest reading will be replaced by a new one. Readings in the memory will not clear away even if power supply is removed.

2. After measurement or when the device stands by, the user can press button M to recall memory. Press button M, the display will show the average value of the latest 3 readings as in Fig.6.

3. Press again, the display will show '01', which means the latest reading, then turns to another screen to show readings as in Fig.7.

4. Press again, the display will show '02', which means the seconds to the latest reading.

Fig. 6

Fig. 7

Memory clearance

When a measurement is finished or when the device is on stand by, hold down button M1 or M2 for at least 5 seconds, the display will show 'Clr' which means the stored reading for M1 or M2 are removed as in Fig.8.

Fig. 8

ERROR AND LOW BATTERY INFORMATION

SYMBOL	POSSIBLE CAUSE	REMEDY
	<p>The cuff is put on wrongly or the tube plug is inserted loosely.</p> <p>Arm/hand movement or talking occurred during measurement</p> <p>The cuff is not inflated to necessary pressure.</p>	<p>Make sure that the cuff is put on correctly and the tube plug is inserted tightly and repeat the measurement.</p> <p>Entirely follow the recommendations in this manual and repeat the measurement.</p>
	<p>The batteries are weak</p>	<p>Replace all 4 batteries with new ones.</p>

MAINTENANCE, STORING, REPAIR AND RECYCLING

1. It is necessary to protect this device against high moisture, direct sunlight, shock, solvent, alcohol and gasoline.
2. Move the batteries if the device is being stored for a long time, and keep the batteries far away from children.
3. Keep the cuff away from sharp objection and don't extend or twist the cuff.
4. This device is not washable. Never immerse the device in water and do not rinse it under the tap. Use only soft and dry cloth to clean the device.
5. Do not serve or maintain the cuff and the device when in use with patient. • Use only soft and dry cloth to clean the device
6. The cuff is sensitive and must be handled with care. You can clean the cuff with damp cloth for daily maintenance. To avoid cross infection when sharing the cuff, you can sterilize the fabric cover of the cuff with tampons moistened by 3% solution of hydrogen dioxide. After long use there will be a partial discoloration on the fabric surface of the cuff. Do not laundry the cuff as well as ironing with a hot flatiron.

WARNING: Under no circumstances may you wash the inner bladder!

7. Since neither the device nor batteries are household waste, follow your local re cycling rules and dispose them at an appropriate collection site.
8. Do not open the device, or delicate electrical components a san intricate air unit could be damaged. If you can not fix the problem using the TROUBLESHOOTING INSTRUCTION, please request service from your dealer.

WARNING: Do not repair the device without manufacturer's authorization. Do no carry out maintenance when using the device.

Caution: Generally, we recommend having the device inspected every 2 years to ensure proper function, accuracy and safety. Please contact your dealer for maintenance.

TROUBLESHOOTING		
SYMPTOM	CHECK POINT	REMEDY
No display after pressing the button 	<p>The batteries have run down.</p> <p>The polarity of installed batteries is wrong.</p> <p>The contact of battery compartment is polluted.</p>	<p>Replace all the batteries with new ones.</p> <p>Install the batteries correctly.</p> <p>Clean the battery terminals with dry cloth.</p>
Inflation stops and re-inflate later.	<p>190mmHg is not enough or the device has moved excessively.</p> <p>Did you talk or move your arm (or hand) during measurement?</p>	<p>See AUTOMATIC INFLATION.</p> <p>Keep quit and silent during the measurement.</p>

SYMPTOM	CHECK POINT	REMEDY
The reading is extremely low or high.	Is the cuff at the same level as the heart. Is the cuff wrapped right? Did you strain your arm during measurement? Did you talk or move your arm (or hand) during measurement?	Make sure that your posture is right. Wrap the cuff correctly. Relax during measurement. Keep quiet and silent during the measurement.
Pulse rate is too low or too high.	Did you talk or move your arm (or hand) during measurement? Did you make measurement right after exercise?	Keep quiet and silent during the measurement. Take measurement again after resting for more than 5 minutes.
The batteries are run down soon.	Faulty batteries are used	Suggest to use alkaline batteries of known manufacturers.

SPECIFICATIONS

Model: INN-006

Size: 110 (L) x 80 (P) X 52 (H) mm

Weight: approximately 215g without batteries

Measuring Method: oscillometric.

Measuring Range:

40 to 180mmHg(DIA,diastolic pressure)

60 to 260mmHg(SYS,systolic pressure)

40 to 160 beats/minute (PUL,pulse rate)

Measuring Accuracy: ± 3 mmHg for static pressure

$\pm 5\%$ of reading for the pulse rate

Inflation: automatic by the pump

Rapid deflation: automatic electronic valve

Batteries: included 4 "AA" x 1,5V

Adapter: optional component, 6V DC 600mA

Memory: 90 memories

Operation Temperature and humidity: +10° C to +40°C; 85% and below (>10%)

Storage/Transportation temperature and humidity: -20°C to +50°C; 85% and below (>10%)

Upper arm circumference: 22~32 cm (standard cuff)

RATED range of CUFF pressure: 300mmHg.

Complete Kit: Main body, standard cuff, 4xAA batteries, gift box, instruction manual, warranty card.

The product is designed with features rated to operate up to an altitude of 2000 m.

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 **1936** The device complies with the European directive 93/42/CEE concerning the medical device and all the other applicable EU norms.

INFORMATION TO USERS according to Legislative Decree No. 49 of March 14, 2014 "2012/19/UE Implementation of the Directive on Waste Electrical and Electronic Equipment (WEEE)" The crossed-out dustbin symbol indicates that the product at the end of its life must be collected separately from other waste. The user should, therefore, take the equipment with the essential components at the end of its useful life to the separate collection center of electronic and electrical waste, or return it to the retailer when purchasing new equipment of equivalent type, in ratio of one to one, or one to zero for devices with larger side less than 25 CM. The separate collection for the decommissioned equipment for recycling, treatment and environmentally compatible disposal contributes to avoid possible negative effects on the environment and human health and promotes recycling of the materials. Improper disposal of the product by the user entails the application of administrative sanctions according to Legislative Decree No. 49 of 14 March 2014.

The batteries used in this device must be disposed of in the special bins at the end of their life. Please inform yourself about the local rules on separate collection of batteries. The correct disposal of batteries helps preventing potentially negative consequences on the environment and human health.


TYPE BF applied parts

IP21


Protection grade


Read carefully the user manual

Temperature range
for storage and
transportation

Air Pressure range
for storage and
transportation

Humidity range for
storage and
transportation

MANUFACTURER'S DECLARATIONS

The model INN-006 Digital Blood Pressure Monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of the model INN-006 Digital Blood Pressure Monitor should assure that is used in such an environment.

EMC GUIDANCE

INN-006 requires special precautions regarding EMC and must be installed and put into service according to the EMC information contained in this user manual; Portable and mobile radiocommunication equipment may affect the operation of INN-006; The appliance or system must not be used near other appliances and, if it is necessary to use it near other appliances, the electro-medical device must be observed to check normal functioning in the configuration in which it is used.

The electro-medical device has been tested and found to comply with the emission and immunity limits of electromedical equipment according to IEC60601-1-2. These limits are designed to provide adequate protection against harmful interference in a typical medical installation. However, there is no guarantee that interference will not occur in a particular installation. If the electro-medical device, interacting with another device, causes or receives detectable interference, the user is invited to limit interference by adopting one or more of the following measures:

1. reorient or relocate the receiving device;
2. increase the distance between the appliances;
3. contact the manufacturer or local technician for assistance.

Guidance and declaration of manufacturer – electromagnetic emissions

INN-006 is intended to be used in the electromagnetic environment specified below. The customer or user of device must assure it's used in this environment.

Emissions tests	Compliance	Electromagnetic environment
Emissions at RF CISPR 11	Group 1	INN-006 uses RF energy only for its internal operation. As a result, RF emissions are very low and are unlikely to cause any interference in nearby electronic devices
CISPR 11 emissions	Classe B	
IEC 61000-3-2 harmonic emissions	CLASS A	
Voltage fluctuations emissions/flicker IEC 61000-3-3	COMPLIES	

INN-006 is suitable for use in all environments, including domestic and those directly connected to a public low-voltage power supply, that supplies buildings used for domestic purposes.

Declaration-Electromagnetic Immunity(Home Healthcare Environment)

Immunity test	IEC 60601 test level	Compliance level
Conducted RF IEC 61000-4-6	3V 150 kHz to 80 MHz 6V in ISM and amateur radio bands between 0.15MHz and 80MHz	3V 150 kHz to 80 MHz 6V in ISM and amateur radio bands between 0.15MHz and 80MHz
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz also meet the requirement of table 9 of 60601-1-2:2014	10 V/m 80 MHz to 2.7 GHz also meet the requirement of table 9 of 60601-1-2:2014
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV,±4 kV,±8 kV,±15 kV air	±8 kV contact ±2 kV,±4 kV,±8 kV,±15 kV air
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines	±2 kV for power supply lines
Surge IEC 61000-4-5	±0.5 kV,± 1 kV line(s) to lines	±0.5 kV,± 1 kV line(s) to lines
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT, 0.5 Cycle at 0°,45°,90°, 135°,180°, 225°, 270°, 315° 0% UT, 1 Cycle and 70% UT, 25/30 cycles sigle phase:at 0° 0% UT,250/300 cycles	0% UT, 0.5 Cycle at 00,45°,90°, 135°,180°, 225°, 270°, 315° 0% UT, 1 Cycle and 70% UT, 25 cycles sigle phase:at 0° 0% UT,250cycles
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m

NOTE: The EUT is the a.c. mains voltage prior to application of the test level.

The following phenomenon is still fulfill the requirement of basic safety and essential performance.

*UT:230V ~/50Hz. The pressure of the EUT is deviation the normal value but the value is still more than 10psi when flow is 4.5l/min.

**UT:230V ~/50Hz. The EUT stop working when adding 0%UT, but the EUT can restore its normal mode automatically.

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30cm (12 inches) to any part of this device, including cables specified by the manufacturer.

Otherwise, degradation of the performance of this equipment could result. Under the test condition specified in immunity, the product can provide the basic safety and essential performance.

If the essential performance is lost or degraded, additional measures are necessary, such as reorienting or relocating the device.

BLOOD PRESSURE RECORDS

DATE	SYSTOLIC (mmHg)	DIASTOLIC (mmHg)	PULSE

innoliving

**Innoliving Spa
Via Merloni, 2/B
60131 Ancona Italy
Tel 071.2133550
www.innoliving.it
MADE IN CHINA**

CE
1936

Rev.05_07.2021