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	<b>JOYTECH Healthcare Co., Ltd.</b> No.365, Wuzhou Road 311100 Hangzhou, Zhejiang Province PEOPLE'S REPUBLIC OF CHINA
	<b>Shanghai International Holding Corp. GmbH (Europe)</b> Eiffestraße 80, 20537 Hamburg, GERMANY
	<b>Innolving Spa</b> Via Merloni, 2/B – 60131 Ancona – Italy www.innolving.it

Thank you for purchasing the INMD-DMT-437 digital thermometer.

This manual has been prepared to ensure safe and efficient use of the DMT-437 thermometer by the end user. The device must be used in accordance with the instructions provided in this manual. It is necessary to read carefully and ensure that you have understood the entire contents of the manual, particularly the sections "IMPORTANT SAFETY INSTRUCTIONS" and "WARNINGS".

## IMPORTANT SAFETY INSTRUCTIONS

### INTENDED USE

This medical device is intended for measuring human body temperature in normal mode orally, rectally, and axillary. The thermometer is a reusable device designed for home or clinical use on individuals of all ages, including children under 8 years old under adult supervision.

### INTENDED USER/INDICATIONS FOR USE

Patients or healthy individuals who wish to measure their body temperature, and healthcare professionals who wish to measure the patient's body temperature. Target patient group: Patients or healthy individuals who need or wish to measure their body temperature.

**CONTRAINDICATIONS:** No contraindications.

### WARNINGS

1. Choking Hazard: the thermometer cap and battery can be fatal if swallowed. Do not allow children to use this device without adult supervision.
2. Do not allow children to use this device without adult supervision.
3. Do not use the thermometer in the ear. Intended use is only for oral, rectal, and axillary readings.
4. Do not place the thermometer battery near extreme heat sources as it may explode.
5. Remove the battery from the device when not in use for an extended period.
6. Using temperature readings for self-diagnosis is dangerous. Consult a physician for interpretation of results. Self-diagnosis can worsen existing health conditions.
7. Do not attempt measurements when the thermometer is wet, as inaccurate readings may occur.
8. Do not bite the thermometer. Doing so may cause breakage and/or injury.
9. Do not attempt to disassemble or repair the thermometer. Doing so may result in inaccurate readings.
10. After each use, disinfect the thermometer, especially if the device is used by multiple people.
11. Do not force the thermometer into the rectum. Discontinue insertion and cease measurement if there is pain. Failure to adhere to this instruction may cause injury.
12. Do not use the thermometer orally after using it rectally.
13. For children aged two years or younger, do not use the device orally.
14. If the unit has been stored at temperatures above 5°C to 40°C (41°F to 104°F), let it reach room temperature (5°C to 40°C) for approximately 15 minutes before use.
15. Avoid using this device near or in contact with other appliances as it may cause improper operation. If such use is necessary, closely observe this device and others to ensure they function properly.
16. RF communication portable equipment (including antenna cables and external antennas) should be used no less than 30 cm (12 inches) from any part of the EM apparatus or EM SYSTEM, including cables specified by the MANUFACTURER. Otherwise, performance degradation of this device may occur.
17. Use in oxygen-rich environments and in the presence of flammable anesthetic mixtures with air, oxygen, or nitrous oxide is not intended.
18. Do not expose the thermometer to direct sunlight or wrap it with cotton, as accuracy will be compromised.
19. EM equipment should not be cleaned and disinfected while in use.
20. Any serious incidents related to the device must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

This digital thermometer provides a quick and highly accurate reading of an individual's body temperature. To better understand its functions and achieve years of reliable results, please read all instructions before use.

This device complies with the following standards:

ISO 80601-2-56 Medical electrical equipment - Part 2-56: Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement.


EN 60601-1-11 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment and meets the requirements of standards.

EN 60601-1-2 (EMC), IEC/EN 60601-1 (Safety). Additionally, the manufacturer is ISO 13485 certified.

## PRODUCT'S DESCRIPTION

Contents:

1. Thermometer
2. User manual
3. Storage case

## PRECAUTIONS

The device's performance may be compromised if one or more of the following conditions occur:

- Operating outside the temperature and humidity range declared by the supplier.
- Storage outside the temperature and humidity range declared by the supplier.

- Mechanical shock (e.g. drop test) or degraded sensor.

- The patient's temperature is lower than the ambient temperature.

Portable and mobile RF communications may affect the device. The device requires special precautions regarding EMC based on the EMC information provided in the attached documents.

## SYMBOLS USED

	Caution		INFORMATION TO USERS according to Legislative Decree No. 49 of March 14, 2014 "2012/19/UE Implementation of the Directive on Waste Electrical and Electronic Equipment (WEEE)". The crossed-out dustbin symbol indicates that the product at the end of its life must be collected separately from other waste. The user should, therefore, take the equipment with the essential components at the end of its useful life to the separate collection center of electronic and electrical waste, or return it to the retailer when purchasing new equipment of equivalent type, in ratio of one to one, or one to zero for devices with larger side less than 25 CM. The separate collection for the decommissioned equipment for recycling, treatment and environmentally compatible disposal contributes to avoid possible negative effects on the environment and human health and promotes recycling of the materials. Improper disposal of the product by the user entails the application of administrative sanctions according to Legislative Decree No. 49 of 14 March 2014.
	Applied part type BF		
	Refer to the instruction for use.		
	The product complies with the requirements of Regulation (EU) 2017/745 MDR on medical devices		
	Manufacturer		
	Authorized Representative in the European Union		
	Importer		
	Manufacturing date		
	Medical device		
	Protection grade		
	Humidity range for transportation and storage		
	Temperature range for transportation and storage		
	Pressure range for transportation and storage		

## TECHNICAL SPECIFICATIONS

Product's description	Flexible digital thermometer	
Model	DMT-437	
Operation mode	Direct (not predictive)	
Measurement range	32.0°C-42.9°C (90.0°F-109.9°F)	
Measurement accuracy	±0.1°C (±0.2°F) at 35.5°C~42.0°C(95.9°F~107.6°F) and for an ambient operating range 18°C~28°C (64.4°F~82.4°F); ±0.2°C (±0.4°F) for other measuring and ambient operating range.	
Display	LCD Display 3 ½ digits	
Memory	Storage of the last measured value	
Battery	1,5V alkaline button cell	
Battery life	Approximately 200 hours of continuous operation or 1 year considering 3 measurements per day	
Dimensions	12.4 cm x 1.9 cm x 1.1 cm (L x W x H)	
Weight	Approx. 11 g including battery	
Expected service life	3 years	
Operating environment	Temperature	5°C ~40°C (41°F ~ 104°F)
	Humidity	15% ~ 95% RH
	Pressure	700hPa ~ 1060hPa
Ambiente di stoccaggio	Temperature	-20°C ~ +55°C (-4°F ~ 131°F)
	Humidity	15% ~ 95% RH
	Pressure	700hPa ~ 1060hPa
Classification	Applied part type BF	
Protection grade	IP27	
Controindications	None	

## UNIT OF MEASUREMENT SELECTION

Temperature readings are available in both Celsius and Fahrenheit scales (°C/°F located in the top right corner of the LCD display). With the unit off, press and hold the ON/OFF button for about 2 seconds to change the current setting

## USE THE DEVICE


1. Press the ON/OFF button located next to the LCD display to turn on the product. The device will emit a short beep, and the screen will display **188.8** followed by the last detected measurement. After displaying the self-test temperature, the thermometer is now in measurement mode.

2. Place the thermometer in the desired measurement site (mouth, rectum, or armpit).

a) Oral use: place the thermometer under the tongue as indicated in the position shown in the figure. Close the mouth and breathe evenly through the nose to prevent the measurement from being influenced by inhaled/exhaled air.

b) Rectal use: Lubricate the silver tip of the probe with petroleum jelly to facilitate insertion. Gently insert the sensor into the rectum for about 1 cm (less than ½ inch).

c) Axillary use: Dry the armpit. Place the probe in the armpit and keep the arm pressed firmly against the side. From a medical standpoint, this method is less accurate for detecting temperature.

3. The degree symbol will flash throughout the measurement. When it stops flashing, an alarm will sound for about 10 seconds, and the reading will appear on the LCD display. It is important to keep the thermometer in place until the audible signal. The measurement will continue even after the signal. For better results in measuring body temperature, it is recommended to keep the probe in the mouth or rectum for about 2 minutes, or under the armpit for about 5 minutes, regardless of the audible signal, and to observe a measurement interval of at least 30 seconds.

Note: The alarm typically sounds faster when the temperature reaches or exceeds 37.8 (100).

4. To prolong battery life, turn off the unit by pressing the ON/OFF button at the end of the measurement. If no buttons are pressed, the unit will automatically turn off after about 10 minutes.

#### TROUBLESHOOTING

Error message	Problem	Solution
<b>L<sub>o</sub></b>	The temperature taken is lower than 32.0°C (90.0°F).	Turn off, wait one minute and take a new temperature via close contact and sufficient rest.
<b>H<sub>i</sub></b>	The temperature taken is higher than 42.9°C (109.9°F).	Turn off, wait one minute and take a new temperature via close contact and sufficient rest.
<b>Err</b>	The system is not functioning properly.	Unload the battery, wait for 1 minute and repower it. If the message reappears, contact the retailer for service.
<b>■</b>	Low Battery: The battery icon flashes, and measurements cannot be taken.	Replace the battery.

#### BATTERY REPLACEMENT

1. Replace the battery when the symbol "■" starts flashing on the LCD display.
2. Remove the battery cover as shown in Figure 1.
3. Gently pull out the compartment containing the battery for about 1 cm (slightly less than 1/2 inch) as illustrated in Figure 2.
4. Use a sharp object such as a pen to remove the old battery. Dispose of the battery in accordance with current regulations. Replace it with a new LR41 1.5V DC button battery. Ensure that the battery is installed with the polarity "+" facing upwards (Figure 3).
5. Slide the compartment with the card and battery inside and close the cover (Figure 4).

Note: The following battery replacement procedure refers to the back of the thermometer.


Figure 1


Figure 2

Figure 3

Figure 4

#### CALIBRATION

The thermometer is initially calibrated at the time of production. If the thermometer is used according to the instructions for use, periodic recalibration is not required. However, we recommend checking the calibration every two years or whenever there is doubt about the clinical accuracy of the thermometer.

If there are signs indicating that the thermometer does not meet the error limits specified in this manual, it is necessary to contact the manufacturer or an authorized technical service.

#### CLEANING AND DISINFECTION

- 1) Immerse the thermometer probe in distilled water for at least 1 minute.
- 2) Use a clean, soft cloth to wipe the thermometer and remove any residue.
- 3) Repeat steps 1 and 2 three times until the probe is visibly clean.
- 4) For complete cleaning and disinfection, use one of the following methods:
  - Method A (High-Level Disinfection): Immerse the thermometer probe in a 0.55% OPA (Ortho-Phthalaldehyde) disinfectant solution for at least 12 minutes at 20°C.
  - Method B (Low-Level Disinfection): Use a clean, soft cloth soaked in 70% medical alcohol and wipe the probe three times, for at least one minute each.
- 5) Repeat steps 1, 2, and 3 to remove OPA residues.
  - Note 1: Rectal use is not recommended for home use as OPA is not readily available outside of hospitals. If rectal measurement is necessary, high-level disinfection is strongly recommended.
  - Note 2: Follow the instructions provided in the OPA manual as a reference.

To prevent damage to the thermometer, please observe the following:

- Do not use benzene, paint thinners, gasoline, or other strong solvents to clean the thermometer.
- Do not attempt to disinfect the thermometer probe (tip) by immersing it in alcohol, OPA, or hot water (over 50°C/122°F) for an extended period.
- Do not use ultrasonic washes to clean the thermometer.

#### RATING LABEL

#### BATTERY INFORMATION

#### EMC GUIDE

The device meets the EMC requirements of the international standard IEC 60601-1-2. These requirements are met under the conditions described in the following table. The device is an electrical medical product and is subject to special EMC precautions that must be published in the user instructions. Portable and mobile HF communication equipment may affect the device. The use of the unit in combination with unapproved accessories may negatively affect the device and alter electromagnetic compatibility. The device should not be used directly adjacent to or between other electrical equipment.

Table 1

EMC Guide – Electromagnetic emissions	
The device is intended for use in the specified electromagnetic environment below. The customer or user of the device must ensure that it is used in such an environment.	
Emissions test	Conformity
Not applicable	Not applicable
CISPR 11 Radiated emission	Group 1 Class B
Harmonic emissions IEC 61000-3-2	Not applicable
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable

Table 2

EMC Guide – Electromagnetic immunity		
The device is intended for use in the specified electromagnetic environment below. The customer or user of the device must ensure that it is used in such an environment.		
Immunity test	IEC 60601 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrostatic transient / burst IEC 61000-4-4	±2 kV for power supply lines 100kHz repetition frequency ±1 kV for input/output lines	NA
Surge IEC 61000-4-5	±0.5 kV, ±1 kV differential mode line-line	NA
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT (100% dip in UT) for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°  0% UT (100% dip in UT) for 1 cycle at 0°  70% UT (30% dip in UT) for 25/30 cycles at 0°  0% UT (100% dip in UT) 250/300 cycles at 0°	NA
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m; 50/60Hz	30 A/m; 50/60Hz
NOTE: U <sub>i</sub> is the AC mains voltage prior to application of the test level.		

Table 3

EMC Guide – Electromagnetic immunity		
The device is intended for use in the specified electromagnetic environment below. The customer or user of the device must ensure that it is used in such an environment.		
Immunity test	IEC 60601 test level	Compliance level
Conducted RF IEC 61000-4-6	3 V <sub>rms</sub> 150 kHz to 80 MHz 6 V <sub>rms</sub> 150 kHz to 80 MHz outside ISMa band	NA
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2,7 GHz	10 V/m
Proximity magnetic fields IEC 61000-4-39	30 kHz, 8 A/m, CW 134, 2 kHz, 65 A/m Pulse modulation 2,1 kHz 13,56 kHz, 7,5 A/m, pulse modulation 50kHz	30 kHz, 8 A/m, CW 134, 2 kHz, 65 A/m pulse modulation 2,1 kHz 13,56 kHz, 7,5 A/m, pulse modulation 50kHz
NOTE 1 At 80 MHz and 800 MHz the higher frequency range applies. NOTA 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. a) The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. the amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz. b) The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,7 GHz are intended to decrease the likelihood that mobile/portable communication equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges. c) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device. d) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.		

Table 4

Recommended separation distances between portable and mobile RF communications equipment and the device				
The device is intended for use in an electromagnetic environment in which radiated therefore disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.				
Frequency MHz	Maximum power W	Distance	IEC 60601 test level	Compliance level
385	1.8	0.3	27	27
450	2	0.3	28	28
710	0.2	0.3	9	9
745				
780				
810				
870	2	0.3	28	28
930	2	0.3	28	28
1720				
28				
1970	2	0.3	28	28
2450	2	0.3	28	28
5240	0.2	0.3	9	9
5500				
5785				
NOTE 1: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.				

**AVVERTENZE**

- This device should not be used in proximity to or directly above other electronic equipment such as cell phones, transmitters, or radio control products. If necessary, the device should be observed to verify normal operation.
- The use of accessories and power cords other than those specified, except for cables sold by the equipment or system manufacturer as spare parts for internal components, may result in increased emissions or decreased immunity of the equipment or system.

	<b>JOYTECH Healthcare Co., Ltd.</b> No.365, Wuzhou Road 311100 Hangzhou, Zhejiang Province PEOPLE'S REPUBLIC OF CHINA
	<b>Shanghai International Holding Corp. GmbH (Europe)</b> Eiffelstraße 80, 20537 Hamburg, GERMANY
	<b>Innoviving Spa</b> Via Merloni, 2/B – 60131 Ancona – Italy www.innoviving.it