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Shenzen Finicare Co., Ltd.

201, No. 50, the 3rd Industrial Park Houting Community,
Shajing Street Bao'an District 518104 Shenzhen
PEOPLE'S REPUBLIC OF CHINA

Riomavix S.L.

Calle de Almansa 55, 1D, 28039 Madrid, SPAIN

Innoliving Spa

Via Merloni, 2/B – 60131 Ancona – Italy
www.innoliving.it

Innoliving Spa

Via Merloni, 2/B – 60131 Ancona – Italy
www.innoliving.it

1. BEFORE USING THE DEVICE**1.1 INTRODUCTION**

Thank you for purchasing the INMD-FC-BP131 blood pressure monitor.

*This unit uses the oscillometric method for blood pressure measurement, detecting the movement of blood through the brachial artery and converting it into a digital reading. This method eliminates the need for a stethoscope during monitoring, making the unit easy to use. The unit automatically stores 2*120 sets of measurement values. You can conveniently review the stored data by pressing the memory button.*

The device includes the following components:

- *Main unit*
- *Arm cuff*
- *User manual*

SYMBOLS USED

	Attention: Indicates a potentially dangerous situation that, if not avoided, could result in death or serious injury.
	Caution: Indicates a potentially dangerous situation that, if not avoided, could cause minor or moderate injury or equipment damage.
	Please read this manual carefully before using the device. Keep it for future reference. Consult your doctor for specific information about your blood pressure.
	Keep dry
	Keep out of direct sunlight
	Fragile, handle with care

 0123	Compliance with EU Regulation (EU) 2017/745 MDR: The product meets the requirements of the Medical Devices Regulation
	Type BF applied part
	Refer to the user manual
SN	Numero di serie
LOT	LOT number
	Manufacturer
EC REP	Authorized Representative in the EU
	Importer
	Date of manufacture
MD	Medical device
UDI	Unique device identifier
	Humidity range for storage and transport
	Temperature range for storage and transport

	Prohibited for use in MRI environments
	<p>INFORMATION TO USERS according to Legislative Decree No. 49 of March 14, 2014 "2012/19/UE Implementation of the Directive on Waste Electrical and Electronic Equipment (WEEE)". The crossed-out dustbin symbol indicates that the product at the end of its life must be collected separately from other waste. The user should, therefore, take the equipment with the essential components at the end of its useful life to the separate collection center of electronic and electrical waste, or return it to the retailer when purchasing new equipment of equivalent type, in ratio of one to one, or one to zero for devices with larger side less than 25 CM. The separate collection for the decommissioned equipment for recycling, treatment and environmentally compatible disposal contributes to avoid possible negative effects on the environment and human health and promotes recycling of the materials. Improper disposal of the product by the user entails the application of administrative sanctions according to Legislative Decree No. 49 of 14 March 2014.</p>

1.2 IMPORTANT SAFETY INSTRUCTIONS

To ensure the proper use of the product, always follow basic safety precautions, including:

- Read all the information in the user manual and any other documentation included in the package before using the device.
- Consult your doctor for specific information about your blood pressure. Self-diagnosis and treatment based on measurement results can be dangerous. Follow your doctor's instructions.
- Use the device only for its intended purpose. Do not use it for other purposes. The device is intended for measuring blood pressure and heart rate in adults. It is not recommended for use on infants either at home or in a medical center.
- Avoid using a mobile phone near the device as it may cause operational malfunction.
- Do not use the device in high-radiation areas to ensure accurate

measurement data.

- Do not disassemble or attempt to repair the device or its components. Avoid using the device in locations with flammable gases (such as an anesthetic gases, oxygen, or hydrogen) or flammable liquids (such as alcohol).
- Avoid using a mobile phone or other devices that emit electromagnetic fields near the unit, as this may cause it to malfunction.
- Note that frequent measurements can cause injury to the PATIENT due to interference with blood flow.
- Remove the batteries if the unit will not be used for three months or more.
- Do not place the CUFF on a wound as it may cause further injury.
- Be cautious of the potential interference with blood flow and possible harm to the PATIENT caused by continuous pressure from the CUFF due to twisting of the connecting tube.
- Adapter requirements: It must meet the following conditions: DC 5V output voltage, 1000mA current, and compliance with IEC 60601-1 and IEC 60601-1-11. Provide two layers of MOPP insulation between the AC input and DC output.
- Avoid installing the batteries with the wrong polarity.
- Replace the batteries with four new ones if they are low.
- If the device is not used for more than three months, remove the batteries to prevent leakage, overheating, breakage, and damage to the body of the blood pressure monitor.
- Intended use: The Arm Blood Pressure Monitor is intended for measuring systolic and diastolic blood pressure and heart rate in adults using a non-invasive technique with an inflatable cuff wrapped around the arm.
- Indication: Intermittent measurement of human blood pressure from the arm.

Contraindications:

- The product should not be used on patients with arrhythmia.
- This product is not suitable for newborns, infants, or women.
- Individuals unable to communicate should not use this product.

Intended users: Professional medical staff and non-professionals.

Patient population: Adults.

Clinical benefits: Measuring blood pressure helps diagnose conditions of the human body.

1.3 WARNINGS AND SAFETY INSTRUCTIONS

- **DO NOT** apply the cuff to an arm being used for intravenous infusion or any other intravascular access, therapy, or arteriovenous (A-V) shunt. Inflating the cuff may temporarily block blood flow, potentially causing harm to the patient.
- The cuff should not be placed on the arm on the side of a mastectomy. For bilateral mastectomy, use the less dominant arm.
- Cuff pressurization may temporarily cause loss of function of monitoring equipment used simultaneously on the same limb.
- Ensure that the use of the **AUTOMATIC** blood pressure monitor does not cause prolonged impairment of blood circulation in the **PATIENT**. Check the affected limb carefully for any signs of reduced circulation.
- If the arm is compressed due to air pressure, it is advisable to loosen the **CUFF** or remove the batteries.

- Avoid touching both the patient and the battery compartment during measurement.
- **WARNING:** Do not use “Luer” connectors. Using “Luer lock” connectors in the tubing may pose a risk of accidental connection to vascular fluid systems, allowing air to be pumped into a blood vessel.
- **WARNING:** Portable RF communication equipment (including peripherals such as antenna cables and external antennas) should be used no less than 30 cm (12 inches) from any part of the blood pressure monitor, including the cables specified by the MANUFACTURER. Otherwise, performance degradation of this equipment may occur.
- **WARNING:** Keep the blood pressure monitor out of the reach of pets, pests, and children.
- **WARNING:** The air tube may pose a strangulation hazard to a child’s or pet’s neck, constituting a life-threatening risk. Do not leave it unattended.
- **WARNING:** Small removable parts (such as batteries, etc.) may be accidentally ingested by a


child, posing a choking or poisoning hazard.

- Keep the product out of the reach of children; use it under adult supervision.
- Report any serious incidents related to the device to the manufacturer and the local competent authority.
- The blood pressure monitor has been clinically investigated according to the requirements of ISO 81060-2:2018+A1:2020


2. OPERATING INSTRUCTIONS

2.1 PRODUCT DESCRIPTION

Main unit


Display


Arm Cuff

Suitable for arm circumference between 22 cm~42 cm (8,7-16,5 inches).

Power supply

- 4 AAA alkaline batteries (not included). Do not use rechargeable batteries.
- Type-C USB cable

USER MANUAL


Please keep the instruction manual for future reference.

2.2 BATTERY INSTALLATION/REMOVAL

2.2.1 Remove the battery compartment cover.

2.2.2 Insert 4 "AAA" batteries, ensuring the positive (+) and negative (-) terminals match the markings in the compartment.

2.2.3 Replace the battery cover.

CAUTION

Ensure you replace the batteries if the low battery symbol appears on the display. Also, avoid leaving low batteries in the compartment, as they might leak and damage the device.

2.3 SETTINGS


2.3.1 Unit settings

Hold the "SET" button for 6 seconds while the unit is off to change the measurement unit (kPa or mmHg). The symbol " $\frac{0}{0}$ ", will appear and blink, indicating the BPM is set to mmHg. Press the "MEM" button to change it to " $\frac{00}{00}$ ", indicating the BPM is set to kPa. Press the "START/STOP" button to confirm the chosen setting.

2.3.2 User settings

Press the "SET" button when the unit is off, and then press the "SET" button again to switch between User 1 and User 2.


2.4 PROPER USE OF THE CUFF

2.4.1 Ensure that the air tube is properly inserted into the main unit.


2.4.2 Remove all clothing from your arm to allow the cuff to fit directly on the skin.

2.4.3 Sit in a chair with your feet flat on the floor. Position your arm on a table to ensure that the cuff is level with your heart.


2.4.4 Thread the arm through the cuff and position the bottom of the cuff approximately 0.8-1.2 inches (2-3 cm) above the elbow.

2.4.5 Securely wrap the cuff around your upper arm using the arm band. Ensure there is a gap of 1-2 fingers between your arm and the cuff. Remove any clothing that might restrict the arm.


2.4.6 Relax your arm, with the palm facing upward and fingers naturally curved. Then, turn on the unit and begin the measurement.

2.5 MEASUREMENT GUIDELINES

2.5.1 Correct Method for Performing the Measurement:

- (1) Ensure your elbow is on the table; sit comfortably with legs uncrossed and back and arm supported.
- (2) Keep the center of the cuff at the same height as the heart.
- (3) Avoid putting clothes inside the cuff.

- (4) Keep palms up and the body relaxed.
 - (5) Sit in a chair with your feet flat on the ground.
- It is recommended to remain relaxed and refrain from talking during the measurement procedure. Allow 5 minutes to elapse before taking the first reading and avoid compressing or restricting the connection tubing.

2.5.2 Ideal Environment for Blood Pressure Measurement:

- (1) Measure in the morning when you are relaxed.
- (2) Ensure that the room temperature is around 20°C.
- (3) Choose a quiet place with minimal surrounding noise.

2.5.3 Measure Blood Pressure at the Same Time Every Day:

Blood pressure fluctuates continuously, and a single measurement may not accurately reflect the situation. It is more reliable to measure blood pressure repeatedly over time. Aim to measure blood pressure daily, preferably in a stable mood, such as early in the morning upon waking up. Try to measure at the same time every day.

Note: Allow at least 2-3 minutes between two consecutive measurements. Depending on individual physiological characteristics, a longer rest period may be necessary.

Ensure that the cuff is around the arm before starting the measurement. Avoid measuring blood pressure while traveling.

Cleaning methods:

- Clean the blood pressure monitor frequently;
- Use a soft, dry cloth to clean the blood pressure monitor; if it is very dirty, you can dampen the cloth with water or neutral detergent, wring it out, and wipe the monitor;
- Disinfect the cuff with cotton wool moistened with 75% alcohol.

2.6 FUNCTIONS

2.6.1. The device features the WHO classification display function.

With the change in blood pressure value, the “” symbol will appear higher for higher blood pressure values (refer to the WHO

- definition of high blood pressure levels - reference table).
- 2.6.2. The device has the function to display the average of the last three measurements (press the "MEM" button to display the value that is the average of the last three measurements).
- 2.6.3. The device has a large display screen.
- 2.6.4. The device allows switching between kPa and mmHg units of measurement and can measure in both modes.
- 2.6.5. The device has a dual memory function and can store 120 sets of measurements for each user, enabling the recall of the normal status of blood pressure.
- 2.6.6. Low battery indicator: under any operating state, the symbol on the display "☐" indicates low battery.
- 2.6.7. Overpressure protection function: when the pressure exceeds 300 mmHg, the device quickly activates automatic power consumption.
- 2.6.8. Automatic shutdown function: if no operation is performed for 2 minutes, the device will automatically shut down.
- 2.6.9. Heartbeat prompting function.
- 2.6.10. Measurement completion prompting function.
- 2.6.11. Error prompting function.

MEASUREMENT range:


Cuff pressure: 0-299 mmHg

Systolic blood pressure (SYS): 45-255 mmHg

Diastolic blood pressure (DIA): 24-222 mmHg

2.7 PERFORM A MEASUREMENT

- 2.7.1. Press the "START/STOP" button. All display symbols will appear on the screen. The cuff will begin to inflate automatically.

2.7.2 Start of measurement. When the pressure reaches the level necessary to stop the blood flow, the monitor stops increasing the pressure and begins to automatically deflate, measuring your blood pressure. The heartbeat symbol flashes when your pulse is detected. When the heart rate symbol flashes, the monitor detects your heart rate and automatically starts calculating your heart rate.

2.7.3 Blood pressure and heart rate are displayed at the end of the measurement. The cuff automatically deflates and all measurement results are saved. The symbol "  " will be displayed if an irregular heartbeat is detected.

2.7.4 Press the "  " button to turn off the unit. The unit will automatically turn off after two minutes if no operation is performed.


NOTE: Inflation or measurement can be stopped at any time by pressing the "  " button.

2.8 USING THE MEMORY FUNCTION

The unit stores blood pressure and heart rate in memory mode after completing each measurement. It can automatically store up to 2x120 sets of measurement values. When more than 2x120 sets are recorded, the oldest record will be automatically deleted to save the latest measured value. The unit also calculates an average reading based on the values of the last 3 measurements.


2.8.1 To read the measurement value

Press the "MEM" button while the unit is off, the unit will enter memory mode and the average value of the user's last 3 measurements will be displayed. Press the "MEM" button again to display the last measurement value, then press the "MEM" button again to check previous measurement values.

2.8.2 To delete memory data

- In memory mode, press and hold the "MEM" button. The unit will display the following symbol to delete all memory records of the current user.
- Press the "START/STOP" button to delete the memory records of the current user. The unit will turn off after deletion.
- Use the memory deletion function with caution.

2.9 BLOOD PRESSURE INFORMATION

2.9.1 When the unit detects an irregular rhythm two or more times during measurement, the irregular heartbeat symbol will appear on the display along with the measured values. An irregular heartbeat is defined as a rhythm that is more than 25% slower or more than 25% faster than the average rhythm detected while the monitor measures systolic and diastolic pressure. If the irregular heartbeat symbol (\sim ) appears in the measurement results, consult your doctor and follow their instructions.


2.9.2 Blood Circulation


Blood circulation is responsible for supplying oxygen to the body.

Blood pressure is the pressure exerted on the arteries. The systolic blood pressure value (higher pressure or upper value) represents the blood pressure produced by the contraction of the heart muscle. The diastolic blood pressure value (lower pressure or lower value) represents the blood pressure produced by the relaxation of the heart muscle.

Systolic Blood Pressure


Diastolic Blood Pressure

2.9.3 Health and Blood Pressure

The incidence of hypertension increases with age. Additionally, lack of exercise, excess body fat, and high levels of cholesterol (LDL) can accumulate inside blood vessels, reducing their elasticity.


Hypertension accelerates atherosclerosis, which can lead to serious conditions such as stroke and myocardial infarction. For these reasons, it is very important to know if blood pressure is within a healthy range. Blood pressure fluctuates from minute to minute throughout the day. Therefore, it is essential to take regular measurements to identify an average blood pressure.

2.9.4 Blood Pressure Classification

After each completed measurement, the display will automatically show the position of the measured pressure on the six segments of the bar indicator, corresponding to the World Health Organization (WHO) blood pressure indicator.

Reference: Journal of Hypertension 1999, Vol 17 No.2

*Note!

When a person's systolic and diastolic pressures fall into different categories, the higher category applies.

2.9.5 Signs of Hypertension

Hypertension can go unnoticed for a long time because it often does not cause obvious symptoms. However, the following conditions can contribute to an abnormal increase in blood pressure:

- Overweight;
- High cholesterol levels;
- Smoking;
- Excessive alcohol consumption;
- Stress and emotional worries;
- Excessive salt intake;
- Lack of physical activity;
- Genetic/hereditary predisposition;
- Presence of preexisting conditions, such as kidney problems or endocrine disorders.

2.9.6 Managing Hypertension

If blood pressure remains consistently high, with values above 140-160 mmHg for systolic and 90-95 mmHg for diastolic, for several days, it is advisable to consult a doctor for a thorough evaluation. You can support the treatment prescribed by your doctor by following these suggestions:

- Lose weight and reduce cholesterol levels
- Limit alcohol consumption
- Reduce salt intake
- Quit smoking
- Maintain a regular physical activity regimen
- Regularly monitor your blood pressure.

3 MAINTENANCE AND CARE

3.1 CLEANING, DISINFECTION AND MAINTENANCE

To keep the blood pressure monitor in optimal condition and protect it from damage, please follow these guidelines:

3.1.1 Clean the blood pressure monitor regularly after each use by the user. Avoid using abrasive or volatile cleaners and never immerse the unit or its components in water.

3.1.2 Use a soft, dry cloth to clean the blood pressure monitor. If it is

heavily soiled, you can dampen the cloth with water or a neutral detergent, wring it out, and clean the monitor.

- 3.1.3 Disinfect the cuff with cotton swabs soaked in 75% alcohol.
- 3.1.4 Store the unit in a safe and dry place. Avoid tightly folding the cuff and protect it from extreme temperatures, humidity, and direct sunlight.
- 3.1.5 Avoid subjecting the unit to strong impacts, such as dropping it on the ground.
- 3.1.6 Remove the batteries if the unit will not be used for three months or longer. Always replace all batteries with new ones at the same time.
- 3.1.7 Use the unit according to the instructions provided in this manual. Use only authorized parts and accessories.
- 3.1.8 Do not attempt to repair or open the device. In case of a defect, contact your local distributor for professional assistance

3.2 CALIBRATION AND SERVICE

The accuracy of this blood pressure monitor has been carefully tested and it is designed for long-term use. As a general recommendation, it is advised to have the unit checked every two years to ensure proper functioning and accuracy. Please consult your local authorized distributor or retailer for assistance.

3.3 ERROR MESSAGES

The following symbols will appear on the display during measurement in case of abnormalities.

Symbol	Cause	Solution
	Anomaly during cuff inflation	Wrap the cuff correctly and tightly, then inflate again.
	Measurement fails	Do not move your arm or body and remain quiet during measurement. Repeat the measurement.
	Low battery	Replace the batteries.

3.4 TROUBLESHOOTING

PROBLEM	CAUSE AND SOLUTION
No power	Replace the batteries
Blank screen	Check that the batteries are installed with the correct polarity.
Measurements are too high or too low	Blood pressure fluctuates constantly. Various factors can influence it, including stress, time of day, and how the cuff is positioned. Refer to the sections "Correct Method for Measurement" and "Taking a Measurement" for more information.

3.5 TECHNICAL DATA

3.5.1 Specifications

Product's description	Upper arm blood pressure monitor	
Model	FC-BP131	
Display	LED Display	
Measurement method	Oscillometric measurement method	
Measurement range	Pressure	0 mmHg – 299 mmHg (0 kPa – 39.9 kPa)
	Heart rate	40 – 180 beats per min.
Measurement Accuracy	Static pressure	±3 mmHg (±3 0.4 kPa)
	Heart rate	Within ±5% of reading
Inflation	Automatic arm cuff inflation	
Deflation	Automatic rapid deflation	
Pressure detection	Semiconductor pressure sensor	
Memory	120 measurements for 2 user groups	
Power supply	4 "AAA" alkaline batteries x 1,5V (not included)	
	USB cable type-C (included)	
	Approximately 300 measurements when using alkaline batteries at the room temperature of 22°C and by using three times a day and inflating to 170 mmHg.	
Automatic shutdown	Within 2 minutes of non-use	
Product weight	Approximately 256.2g (excluding batteries)	
Dimensions	110mm x 140mm x 67mm	
Cuff circumference	22 ~ 42 cm	
Operating environment	Temperature	10°C ~ 40°C (50°F-104°F)
	Humidity	15% ~ 85% RH
Transport and Storage environment	Temperature	-20°C ~ +55°C (50°F-140°F)
	Humidity	10%-95% RH
	Atmospheric pressure	70-106 kPa
Classification	Internal power supply appliance. Applied part (cuff) type BF 	
Accessories	Rigid cuff, user manual and USB cable	
Service life	5 years	

o enhance performance, these specifications may be subject to change without notice. The device, accessories, and packaging must be disposed of properly at the end of use to reduce risk to the patient or user to an acceptable level.

3.5.2 IEC 60601-1-2:2014/AMD1:2020 EMC and EMF Equipment and Systems Identification, Marking and Product Documentation for Class B Equipment.

USER INSTRUCTIONS:

The EMC or EMF equipment/system is suitable for healthcare use in home and similar environments.

WARNING: Avoid using the device near active high-frequency surgical equipment and MRI shielded rooms, where electromagnetic disturbances are high.

WARNING: Do not place the device near other equipment or stack it with them, as this may cause malfunctions. If necessary, monitor the operation of both devices to ensure they operate correctly.

Warning: The use of accessories, transducers, and cables other than those supplied by the manufacturer may increase electromagnetic emissions or reduce electromagnetic immunity, leading to operational issues.

WARNING: Portable RF communication equipment, such as antenna cables and external antennas, should be used at a distance of at least 30 cm (12 inches) from any part of the Electronic Arm Blood Pressure Monitor (Model: FC-BP131). A shorter distance may compromise the device's performance.

IF APPLICABLE, provide a list of all cables, with maximum length, transducers, and other replaceable ACCESSORIES provided by the RESPONSIBLE ORGANIZATION that may affect the compliance of ELECTROMEDICAL EQUIPMENT or SYSTEMS with the requirements of Clauses 7 (EMISSIONS) and 8 (IMMUNITY). ACCESSORIES may be specified generically (e.g., shielded cable, load impedance) or specifically (e.g., by MANUFACTURER and TYPE REFERENCE).

IF APPLICABLE, describe the performance of the ELECTROMEDICAL EQUIPMENT or SYSTEMS considered ESSENTIAL and explain what the OPERATOR can expect in case of loss or degradation of such performance due to ELECTROMAGNETIC DISTURBANCES (the defined term “ESSENTIAL PERFORMANCE” is not required).

RATING LABEL

3.6 EMC GUIDE

1. All necessary instructions to ensure the maintenance of BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic interference for the entire expected service life.
2. Manufacturer's indications and statements regarding electromagnetic emissions and immunity.

Table 1

Guidance and manufacturer's declaration – electromagnetic emissions	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations/flicker emissions IEC 61000-3-3	Compliant

Table 2

Guidance and manufacturer's declaration – electromagnetic immunity		
Immunity test	IEC 60601-1-2 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	Power supply lines: ±2 kV, 100kHz repetition frequency	Power supply lines: ±2 kV, 100kHz repetition frequency
Surge IEC 61000-4-5	Line(s) to line(s): ±1 kV	Line(s) to line(s): ±1 kV
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT 1 cycle and 70% UT, per 25/30 cycles single phase at 0° 0% UT 300 cycle	0% UT 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT 1 cycle and 70% UT, per 25/30 cycles single phase at 0° 0% UT 300 cycle
Power frequency magnetic field IEC 61000-4-8	30 A/m; 50/60Hz	30 A/m; 50/60Hz
Radiated RF IEC 61000-4-3	10 V/m Da 80 MHz a 2.7 GHz 80% AM a 1kHz	10 V/m Da 80 MHz a 2.7 GHz 80% AM a 1kHz
Proximity magnetic fields IEC 61000-4-39	30 kHz: 8 A/m 134.2 kHz: 65 A/m 13.56 MHz: 7.5 A/m	30 kHz: 8 A/m 134.2 kHz: 65 A/m 13.56 MHz: 7.5 A/m
NOTE: UT is the AC mains voltage prior to application of the test level.		

Table 3

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Table 4

Guidance and manufacturer's declaration – electromagnetic immunity		
Test frequency	Modulation	IMMUNITY TEST LEVEL (A/m)
30 kHz	CW	8
134,2 kHz	Pulse modulation a 2,1 kHz	65 ^b
13,56 MHz	Pulse modulation a 50 kHz	7,5 ^b
a) The carrier shall be modulated using a 50% duty cycle square wave signal.		
b) r.m.s., before modulation is applied.		

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference.
- (2) This device must accept any interference received, including interference that may cause undesired operation.

NOTE: This device has been tested and complies with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in residential installations. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment causes interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Unauthorized modifications to this device could void the user's authority to operate it.

4 WARRANTY

The warranty for this blood pressure monitor is 3 years from the date of purchase. Warranty obligations are outlined in the warranty card for the buyer.


WARNING:

- Do not modify this device without manufacturer authorization.
- All major device maintenance operations must be performed by an authorized service center or distributor. There are no user-repairable parts inside; contact an authorized representative or the manufacturer

Shenzen Finicare Co., Ltd.

201, No. 50, the 3rd Industrial Park Houting Community,
Shajing Street Bao'an District 518104 Shenzen
PEOPLE'S REPUBLIC OF CHINA

Riomavix S.L.

Calle de Almansa 55, 1D, 28039 Madrid, SPAIN

Innoliving Spa

Via Merloni, 2/B – 60131 Ancona – Italy
www.innoliving.it

Innoliving Spa

Via Merloni, 2/B – 60131 Ancona – Italy
www.innoliving.it

innoliving
easy tech. easy life

Innoliving Spa
Via Merloni, 2/B
60131 Ancona Italy
Tel 071.2133550
www.innoliving.it

MADE IN CHINA

CE
0123

Rev.00_06.2024