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FLEXIBLE DIGITAL THERMOMETER

INSTRUCTION MANUAL

Thank you for purchasing the Innolving digital thermometer INN-018. Before use, it is advisable to carefully read the following instructions and keep this manual safe for future references. This digital electric thermometer is suitable to take body temperature for clinical and hospital use as well as for home use. Operator shall be at least 11 years old and patient can be operator. Patient's weight must be over 2.5kg.

IMPORTANT SAFETY INSTRUCTIONS

- Do not allow the device to come into contact with hot water.
- Do not expose the device to high temperatures or direct sunlight.
- Do not drop the thermometer, this product is not shock resistant.
- Do not attempt to tamper with or repair the product by yourself.
- Do not use thinners, petrol or benzene to clean the device.
- Clean only with water or disinfectant.
- The thermometer contains small parts (such as the battery, and the battery compartment) that can be swallowed by children. Never leave the product unattended.
- Do not bend the thermometer tip.
- If the ambient temperature is above 35 °C or 95 °F, immerse the tip in cold water for about 5-10 seconds before measuring body temperature.
- In case of persistent fever, especially in children, consult a doctor.
- Do not use the thermometer near strong electromagnetic fields; keep it away from any radio system or mobile phone.
- Do not immerse the thermometer more than 15cm in water and do not leave it submerged for more than 30 minutes.

HOW TO USE

Press the **ON/OFF** button to turn on the device, the appliance will emit a short acoustic signal indicating that the thermometer is operational. At the same time, the thermometer runs to self-check test, during which all the digital segments appear on the LCD. Then the last measurement is shown and then the display will show the symbol "L" in the middle and the "°C" in the upper right corner. The device is ready to measure the temperature when the "°C" symbol will start flashing. If the ambient temperature is below 32°C or 89.6°F, the display will show "L °C" or "L °F". If the temperature is higher than 42 °C or 107.6 °F, the display will show "H °C" or "H °F". During the reading phase, the display shows the temperature variation and the "°C" or "°F" symbol flashes in the top right corner. The measurement is completed when the device reaches a constant temperature. The temperature is considered constant when the oscillations are lower than 0.1 °C for about 16 seconds. As soon as the appliance detects the constant temperature, the thermometer emits 4 beeps, the "°C" or "°F" symbol stops flashing and the display shows the detected temperature. At the end of the measurement, turn off the thermometer by keeping the ON/OFF button pressed. If not, the thermometer will automatically shut off after 10 minutes of inactivity.

MEMORY FUNCTION

The thermometer stores in memory the last detected temperature. The measured temperature is automatically displayed for about two seconds after switching on of the device immediately after the self-check test. If a new measurement is done this will be stored replacing the previous one and will be displayed on the next switch on.

MEASUREMENT METHODS

It is important to remember that the body temperature reading depends on the measurement method. For this reason, it should always be taken into consideration how the temperature is detected in order to have a correct reading of the values:

Rectal Use: This is the most accurate method from a medical point of view, because it is closer to the actual body temperature.

Insert the probe into the rectum (2cm max.) And wait for the beep. The measurement time is approximately 40-60 seconds.

Axillary Use: In this way, you will have a temperature measurement that can change from 0.5 °C to 1.5 °C from the actual body temperature. The measurement time is approximately 80-120 seconds.

Oral Use: The oral temperature ranges from 0.3 °C to 0.8 °C compared to the rectal temperature. To ensure that the reading is as accurate as possible, place the tip of the thermometer under the tongue and keep the mouth tightly closed and breathe through the nose. Do not eat or drink anything before the measurement. The measurement time is approximately 50 to 70 seconds.

Note: the rectal mode is considered the most accurate method to identify the body temperature and it is recommended to prolong the measurement time of 3 minutes after the acoustic signal.

BATTERY REPLACEMENT

When the symbol **■** or **□** flashes on the display, it means that the battery level is low and needs to be replaced. Remove the battery compartment and replace the battery with one of the same type, the "+" sign up and "-" sign down. Remove the battery from the device, in case it is going to be unused for a long period of time.

CLEANING AND DISINFECTION

The best way to clean the thermometer tip is applying a disinfectant (e.g. 70% medical alcohol) on a damp cloth and then proceed to clean and disinfect the probe. It is advisable to disinfect the probe before each use. Remove the battery from the device, in case it is going to be unused for a long period of time.

TECHNICAL SPECIFICATION

Measuring Range: (32.0~ 42.0)°C/(89.6~107.6) °F

Measuring Accuracy: +/- 0.1°C/0.2°F (35.5°C~42.0°C/95.9°F~107.6°F) +/-0.2 °C/0.4 °F (32.0°C~35.5°C/89.6°F~95.9°F)

Storage Conditions: (-25~55)°C, ≤95%RH

Ambient Temperature during use: (5~35)°C, ≤80%RH

Minimum Detection: 0.1°C/0.1°F

Battery: Alkaline battery LR41, 1.5 V, minimum duration 100 hours in continuous operation.


Weight: about 19g

Atmospheric Pressure: 700~1060hPa

SYMBOLISM:

- **■** or **□**: Low battery
- "L °C" or "L °F": The temperature is lower than 32 °C or 89.6°F
- "H °C" or "H °F": The temperature is higher than 42 °C or 107.6°F

PRODUCT DESCRIPTION

Legal requirements and guidelines

This product complies with the European directive on medical devices 93/42 / EEC and has the CE mark. The device also complies with the specifications of the European standard EN 12470-3 Clinical thermometers - Part 3: performance of compact electric thermometers (non-predictive and predictive) with the maximum device. The CE mark confirms that it is a medical device with a measuring function under the Medical Device Act which has been subject to a conformity assessment procedure. A Notified Body confirms that this product complies with all applicable laws.

Calibration check

This thermometer is calibrated during production. If the product is used according to the instructions in this manual, a periodic calibration is not necessary. If there are signs indicating that the thermometer does not respect the error limits indicated in this manual, please contact the manufacturer or the authorized technical assistance service.

EMC Guide

INN-018 needs special precautions regarding EMC and must be installed and commissioned in accordance with the EMC information contained in this manual. Portable and mobile radiocommunication devices can affect the operation of INN-018. The device or the system must not be used near to other devices; if it were necessary to use near to other devices, the medical device must be supervised to check normal operation under the setting in use. The medical device has been tested and found compliant with the emission and immunity limits of medical devices according to IEC60601-1-2:2014. These limits are fixed to provide proper protection against harmful interference in a typical medical installation. However, there's not any guarantee that interference will not occur in a particular installation. If the medical device, interacting with another device, causes or receives detectable interferences, the user is invited to limit the interference by adopting one or more of the following measures:

1. to redirect or to reposition the receiving device
2. to increase the distance between the devices
3. to contact the manufacturer or local technician for assistance

Guide and declaration of the manufacturer - electromagnetic emissions

INN-018 is designed to work in the electromagnetic environment specified below. The customer or the user of the appliance should ensure that it is used in this environment

Emission tests	Compliance	Electromagnetic environment
Emissions at RF CISPR 11	Group 1	INN-018 uses RF energy only for its internal operation. As a result, RF emissions are very low and is likely to cause no interference in nearby electronic devices.
CISPR 11 Emissions	Class B	INN-018 is suitable for use in all environments, including domestic environments and those directly connected to a public low-voltage power supply that supplies buildings used for domestic purposes
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

Guide and declaration of the manufacturer - electromagnetic immunity

INN-018 is designed to work in the electromagnetic environment specified below. The customer or the user of the appliance should ensure that it is used in this environment

Immunity tests	Test level of IEC 60601	Conformity level	Electromagnetic environment
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 8kV contact discharge +/- 2; 4; 8; 15 kV air discharge	+/- 8kV contact discharge +/- 2; 4; 8; 15 kV air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
High-frequency electromagnetic field (50/60Hz) IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should have characteristic levels of a typical location in a commercial or hospital environment.

Guide and declaration of the manufacturer - electromagnetic immunity

INN-018 is designed to work in the electromagnetic environment specified below. The customer or the user of the appliance should ensure that it is used in this environment. Portable and mobile RF communications equipment should not be used closer to any part of INN-018, including cables, than the recommended separation distance calculated with the equation applicable to the transmitter frequency.

Immunity test	Test level of IEC 60601	Conformity level	Recommended separation distance d:
Irradiated RF IEC 61000-4-3	3 V/m and 10 V/m from 80 MHz to 2,7 GHz	3 V/m and 10 V/m	d= 30 cm
Immunity to the electromagnetic fields generated by other communication devices RF wireless IEC 61000-4-3	TETRA 400 380 - 390 MHz	27 V/m 27 V/m	d= 30 cm
	GMRS 460 FRS 460 430 - 170 MHz	28 V/m 28 V/m	
	LTE Band 13, 17 704 - 787 MHz	9 V/m 9 V/m	
	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5 800 960 MHz	28 V/m 28 V/m	
	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 5 1700 - 1990 MHz	28 V/m 28 V/m	
	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 70 2400 - 2570 MHz	28 V/m 28 V/m	
WLAN 802.11 a/n 5100 - 5800 MHz	9 V/m 9 V/m		


RATING LABEL:

Application parts type BF.

Read carefully the instruction manual.

Read the user manual for important for important cautionary information, such as warnings and safety norms.

The product complies with the EU directive 93/42/CEE regarding the medical devices, then with all the applicable EU norms.

Lot number (first 2 digits shows the week of production, the last 2 shows the year of production)

The batteries used in this device must be disposed of in the special bins at the end of their life. Please inform yourself about the local rules on separate collection of batteries. The correct disposal of batteries helps preventing potentially negative consequences on the environment.

IP22 Protection against the penetration of solid foreign objects larger than ϕ 12.5 mm and against vertically falling drops when the the enclosure is tilted at any angle up to 15°.

 INFORMATION TO USERS according to Legislative Decree No. 49 of March 14, 2014 "2012/19/UE Implementation of the Directive on Waste Electrical and Electronic Equipment (WEEE)". The crossed-out dustbin symbol indicates that the product at the end of its life must be collected separately from other waste. The user should, therefore, take the equipment with the essential components at the end of its useful life to the separate collection center of electronic and electrical waste, or return it to the retailer when purchasing new equipment of equivalent type, in ratio of one to one, or one to zero for devices with larger side less than 25 CM. The separate collection for the decommissioned equipment for recycling, treatment and environmentally compatible disposal contributes to avoid possible negative effects on the environment and human health and promotes recycling of the materials. Improper disposal of the product by the user entails the application of administrative sanctions according to Legislative Decree No. 49 of 14 March 2014.

Involving Spa
Via Merloni, 9 - 60131 Ancona Italy
Tel: 071 2133550 - www.involving.it
MADE IN CHINA

Rev.01_11.2018