innoliving

easy tech, easy life

**BILANCIA DIAGNOSTICA
BODY FAT & BODY ANALYZER**

INN-117
MANUALE D'USO
USER MANUAL

BILANCIA DIAGNOSTICA BODY FAT & BODY ANALYZER

I

MANUALE D'ISTRUZIONE

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THE DIAGNOSTIC SCALE BODY FAT & BODY ANALYZER

INSTRUCTION MANUAL

Thank you for choosing the diagnostic scale Body Fat & Body Analyzer INN-117. In order to ensure the correct usage of the scale, it is advisable to read carefully this instruction manual. It is also advisable to save this manual for a future use.

WARNINGS

- In order to ensure the correct usage of the scale, it is advisable to read carefully this instruction manual. An incorrect use can damage the device.
- Use this product only for its intended purpose, as described in this instruction manual such as a bathroom scale for domestic use only. Any other use should be considered improper and dangerous.
- Packaging components (plastic bags, cardboard, polystyrene etc.) must be kept out of the reach of children as they have the potential to be dangerous and their disposal should be carried out according to the current regulations.
- This scale should not to be used to weigh objects or substances in commercial transactions, to produce medicines, to calculate tolls, tariffs, taxes, premiums, fines, remunerations, indemnities or fees of a similar kind determined by weight.
- This scale is to be used for keeping one's own personal weight under control, and is not to be used for diagnostic purposes or for medical treatment. Values which may seem unusual or outside the norm should always be discussed with your doctor. In order to avoid danger to the user of this device, the doses of any drug or therapy prescribed by one's own doctor should in no circumstances be modified on the basis of data calculated by this device.
- Avoid overloading the scale, it may be damaged.
- Do not open or tamper with the device: repairs must be performed only by authorized personnel; otherwise the warranty expires.
- Keep the scale platform dry to avoid slipping.
- Do not place the balance on a soft surface (like a carpet) otherwise the weight measurement may be inappropriate.
- Do not place any objects on the scale when not in use.
- Do not jump on the scale.
- It is prohibited the use of the product to people with medical implants (e.g. pacemakers), as the scale may interfere with the medical device.
- Do not use if you are pregnant.
- Caution! Do not step on the scale with wet feet.
- Keep batteries and scale out of reach of children.
- The detection of parameters of children under the age of 10 years may not be correct.
- The scale can be used only for the detection of weight, without disclosing the percentage of water, body fat, muscle mass and other parameters in children of any age. The minimum weight limit is 2kg.
- Unreliable results in measurement of water, body fat or muscle mass could occur in case of professional athletes and people at risk of edema.
- We reserve the right to make changes in the design of the scale, without affecting the safety, performance or functionality of the device.

TECHNICAL SPECIFICATIONS

- Capacity: 150kg
- Division: 100g
- Unit of measure: Kg/lb/st
- Display of body weight, body fat, body water, and muscle ratio.
- Big backlighted display
- Low battery indicator
- Overload indicator
- Power: 2 AAA batteries included.
- Work temperature 5 °C ~ 40 °C

PRINCIPLE OF OPERATION

The principle behind the operation of this diagnostic scale is bioelectrical impedance analysis (BIA). The measurement of the percentage of the different fabrics occurs in a few seconds by means of an electric current, not perceptible to the human being, safe and harmless.

It is important that during the measurement that you consider as follows: measuring the percentage of body fat should only be performed in bare feet; it is advisable to slightly moisturize the soles of the feet. Soles completely dry could be the cause of an unsatisfactory result, because of the low electrical conductivity.

Remain fixed during the measurement.

Make sure that the scale is placed on a flat, stable surface. Do not place the balance on a carpet or on soft surfaces.

To ensure a correct measurement, it is necessary to weigh yourself without clothing or shoes, and before eating.

To ensure accurate measurement of weight, you should not weigh yourself before two hours after waking.

Step on the scale. Do not lean on anything, stay straight and not move.

TURN THE SCALE ON

Insert the AAA batteries into the battery compartment according to polarity. Close the compartment. The scale is ready to operate.

Select the unit of measurement using the switch.

FUNCTION WEIGHT ONLY

Step on the scale, the scale will automatically turn on.

After a few seconds, the scale will detect the weight.

AUTOMATIC POWER OFF

The scale switches off automatically when not in use after 20 seconds .

After detecting the weight, the scale switches off automatically after 6 seconds.

OVERLOADING

If the scale is overloaded, the display will show "EEEE".

BATTERIES REPLACEMENT

If the display shows "Lo" replace the batteries.

DISPLAY SYMBOLS

Display	Significato
	Man
	Woman
CM (FT;in)	Height
Age	Age
FAT	Fat
TBW MUS	Water muscle ratio
BONE %	Bone weight percentage

DATA FOR THE DETECTION OF FAT, WATER, MUSCLE, BONE

You must insert your personal data before making the measurement of various parameters related to the physical form (fat, water, muscle, bone).

Press the SET button, the first of 9 users (P0) will flash on the display. If you want to enter the data in this user, press SET, or select another profile using the arrows and confirm with SET.

Once you confirmed the user, press SET, a flashing  man or  woman will appear, select using the arrows and confirm with SET.

Now the age will flash, insert your age using the arrows and confirm with SET.

Then, the height in cm will flash on the display, set the correct value and confirm with SET. After confirming this last parameter (height, cm), the display will show 0.0 KG and the user name, for example, P0. Start the measurement stepping on the scale. Stand still during the measurement until the data are not detected.

Warning: during the measurement there should not be any contact between the two feet, calves, shins or thighs. Otherwise, the measurement can not be performed properly. Initially, the display will show your weight, then after few seconds the following information will be shown one after the other:

Display	Meaning
FAT	Fat
TBW	Water
MUS	Muscle Ratio
BONE	Bone weight

The scale gives you a general indication of your fitness state according to the following parameters

Age	Woman		Man		RESULTS	
	Fat	Water	Fat	Water	Indication	Symbol
≤30	4.0-20.5	66.0-54.7	4.0-15.5	66.0-58.1	Underweight	-
	20.6-25.0	54.6-57.6	15.6-20.0	58.0-55.0	Normal	0
	25.1-30.5	51.5-47.8	20.1-24.5	54.9-51.9	Overweight	+
	30.6-50.0	47.7-34.4	24.6-50.0	51.8-34.4	Obese	++
≥30	4.0-25.0	66.0-51.6	4.0-19.5	66.0-55.4	Underweight	-
	25.1-30.00	51.5-48.1	19.6-24.0	55.3-52.3	Normal	0
	30.1-35.0	48.0-44.7	24.1-28.5	52.2-49.2	Overweight	+
	30.1-50.0	44.6-34.4	28.6-50.0	49.1-34.4	Obese	++

OTHER SYMBOLS ON THE DISPLAY

WARNING: If you get an error message on the display, wait until the device turns off. Turn on the scale again to perform the measurement again. If the error message on the display doesn't disappear, remove the batteries from the compartment, after 10 seconds insert them again and turn on the device.

Symbols	Meaning
8888	User selection and personal data.
oooo	Data analysis after measurement
----	-Improper placement of the scale. Place the scale in a flat and stable surface. - You are not on the scale with bare feet. Remove your clothes and shoes before stepping on the scale.
LO	Low battery indicator. Replace the batteries.
ErrH	High % of body fat indicator. The percentage of body fat is too high. It is advisable to adopt a low-calorie diet, do exercise and pay more attention to your health.
Errl	Low % of body fat indicator. The percentage of body fat is too low. It is advisable to adopt a more nourishing diet and pay more attention to your health.
EEEE	Overload indicator. The weight exceeds the maximum threshold of measurement. Step off the scale to avoid damaging the product.

In case of malfunctioning:

- Check that the batteries are inserted correctly
- Check that the scale is on a flat, stable surface

MAINTENANCE

Wipe with a damp cloth: avoid any infiltration of water. Never use chemical detergents. Do not immerse in water. Do not step on the scale with wet feet. Do not store the scale in an upright position as this facilitates the batteries run out. As a precision device, the balance must always be treated with care. Do not drop it and do not jump over it. Do not open or tamper with the device: repairs must be performed only by authorized personnel; otherwise the warranty expires. If the product is not used for a long period of time, it is advisable to remove the battery, as a possible leakage could damage your equipment.

Also because the quality of the batteries vary from model to model, the product warranty does not cover damage caused by a possible leak. When the scale is not in use, do not leave weights resting on the platform. Keep it away from dust and moisture; do not expose to the sun or to excessive temperatures.

RATING LABEL:

innoliving

Innoliving Spa
Via Merloni,9 60131 Ancona Italy

INN-117 DIGITAL BATHROOM SCALE
MAX: 150 KG DIV:100 G

3V  (2XAAA batt.)

HOUSEHOLD USE ONLY

MADE IN CHINA
SOLO USO DOMESTICO

LOT

Rev. 00_10.2015

This product complies with all applicable European directives

INFORMATION TO USERS according to Legislative Decree No. 49 of March 14, 2014 " 2012/19/UE Implementation of the Directive on Waste Electrical and Electronic Equipment (WEEE) "The crossed-out dustbin symbol indicates that the product at the end of its life must be collected separately from other waste. The user should, therefore, take the equipment with the essential components at the end of its useful life to the separate collection center of electronic and electrical waste, or return it to the retailer when purchasing new equipment of equivalent type, in ratio of one to one, or one to zero for devices with larger side less than 25 CM. The separate collection for the decommissioned equipment for recycling, treatment and environmentally compatible disposal contributes to avoid possible negative effects on the environment and human health and promotes recycling of the materials. Improper disposal of the product by the user entails the application of administrative sanctions according to Legislative Decree No. 49 of 14 March 2014

The batteries used in this device must be disposed of in the special bins at the end of their life. Please inform yourself about the local rules on separate collection of batteries. The correct disposal of batteries helps preventing potentially negative consequences on the environment.

innoliving
easy tech, easy life

Innoliving Spa
Via Merloni, 9
60131 Ancona Italy
Tel 199.188.335
www.innoliving.it

MADE IN CHINA

CE

Rev.00_10.2015