innoliving

easy tech, easy life

**BILANCIA DIGITALE
BODY FAT & BODY WATER**

INN-140
MANUALE D'USO
USER MANUAL

MANUALE D'ISTRUZIONE

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INSTRUCTION MANUAL

Thank you for purchasing Body Fat & Body Water digital scale INN-140.

Before use, read this manual carefully and store it for future use.

WARNINGS

- For a proper use of the product, please read the following instruction manual. Incorrect use can cause damage to the device.
- Use this product only for its intended purpose, as described in this instruction manual such as a bathroom digital scale for household use only. Any other use should be considered improper and dangerous.
- Keep all packaging materials (plastic bags, cardboard, polystyrene etc.) out of the reach of children as they are potentially dangerous and must be disposed of as required by law.
- This scale should not to be used to weigh objects or substances in commercial transactions, to produce medicines, to calculate tolls, tariffs, taxes, premiums, fines, remunerations, indemnities or fees of a similar kind determined by weight.
- This scale is to be used for keeping one's own personal weight under control and is not to be used for diagnostic purposes or for medical treatment. Values which may seem unusual or outside the norm should always be discussed with your doctor. In order to avoid danger to the user of this device, the doses of any drug or therapy prescribed by one's own doctor should in no circumstances be modified on the basis of data calculated by this device.
- Avoid overloading the scale, it may be damaged.
- Do not open or tamper with the device: repairs must be performed only by authorized personnel; otherwise the warranty expires.
- Keep the scale platform dry to avoid slipping.

- Do not place the balance on a soft surface (like a carpet) otherwise the weight measurement may be inappropriate.
- Do not place any objects on the scale when not in use.
- Do not jump on the scale.
- Do not overload the scale.
- Warning! Do not step on the scale with wet feet.
- Keep batteries and scale out of reach of children.
- The detection of parameters in children under the age of 10 years may not be correct.
- The balance can be used only for the detection of weight, without disclosing the percentage of water, body fat, muscle mass and other parameters in children of any age. The minimum weight is 2.5kg
- The manufacturer reserves the right to make changes on the design of the scale, without affecting the safety, performance or functionality of the device.
- The detection of body fat is not indicated in the following cases:
 - Children under the age of 10 and adults over 85 years.
 - Persons with fever, oedema or osteoporosis.
 - People on dialysis treatment.
 - People with pacemakers.
 - People taking cardiovascular medicament.
 - Pregnant women
 - Athletes who perform more than 5 hours of high-intensity sports activities every day
 - People with heart rate below 60 beats / minute
- Avoid contact with metal objects (rings, chains), there is a risk of short circuit. As a result of the short circuit, the battery can heat up considerably or even catch fire causing burns.

- For safety reasons, the batteries must be covered with adhesive strips during transport. If the battery fluid leaks out and gets in contact with eyes, rinse immediately and seek medical advice.

TECHNICAL SPECIFICATIONS

- Capacity: 150kg
- Division: 100g
- Unit of measure: Kg/lb
- Number of users: 10
- Measure of weight, fat mass and percentage of water
- Body fat measuring unit: 0.1%
- Body fat measuring range: 4.0%-50.0%
- Body water measuring range 35.0%-85.0%
- Wide backlighted display
- Low battery indicator
- Overload indicator
- Power: 2 batteries CR2032 included.
- Operating temperature 5 °C ~ 40 °C

PRINCIPLE OF OPERATION

The principle behind the operation of this diagnostic scale is bioelectrical impedance analysis (BIA). The measurement of the percentage of the different tissues appears in a few seconds by means of an electric current, not perceptible to the human being, safe and harmless.

During the measurement, it is important to consider as follows: measuring the percentage of body fat should only be performed with bare feet; it is advisable to slightly moisturize the soles of the feet. Soles completely dry could be the cause of not satisfactory result, because of the low electrical conductivity. Remain fixed during the measurement.

Make sure that the scale is placed on a flat, stable surface. Do not place the scale on a carpet or on soft surfaces.

To ensure correct measurement, it is necessary to weigh yourself without clothes or shoes, and before eating.

To ensure an accurate measurement of weight, you should not weigh yourself

OPERATING INSTRUCTIONS

Insert the included batteries CR2032 into the compartment according to polarities.

Close the compartment. The scale is ready to operate. Select the unit of measure using the switch at the rear of the balance.

WEIGHTING FUNCTION

Step on the scale, the device will turn on automatically.

After a few seconds your weight will be displayed.

AUTOMATIC SWITCHING OFF

Once measured the weight, the scale will switch off automatically after a few seconds.

OVERLOADING

If the scale is overloaded, the display shows "EEEE"

BATTERY SUBSTITUTION

If the display shows "Lo", replace batteries.

If the scale is overloaded, the display will show "EEEE".

BATTERIES REPLACEMENT

If the display shows "Lo" replace

DISPLAY SYMBOLS

Display	Significato
	Man
	Woman
CM (FT;In)	Height
Age	From 10 to 80
FAT	Fat
TBW	Water
%	Percentage

DATA FOR THE DETECTION OF FAT, WATER, MUSCLE, BONE

You must enter your personal details before making the measurement of various parameters related to the physical form (fat and water).

Press the SET button, the first of 10 users (P1-02-03-04- ...) flashes on the display. If you want to enter data in this user, confirm with SET, otherwise choose another profile using the arrow keys and confirm with SET.

After confirming the user to use, by pressing the SET, a  man or  woman will flash, select using the arrows and confirm with SET.

Now the Height in cm flashes, enter your height using the arrows and confirm with SET.

Then the age will flash, set the correct value and confirm with SET.

After confirming the last parameter (age), the display will show 0.0 KG and the user name, for example, P1. To start the measurement, step on the scale. Remain fixed during the measurement until the data is not detected.

Warning: during the measurement there should not be any contact between the two feet, calves, shins or thighs. Otherwise, the measurement cannot be performed properly.

Initially, the display will show your weight, then after a few seconds the following information will be shown one after the other

Display	Meaning
FAT	Fat
TBW	Water

This scale gives you a general indication of your fitness state with the following parameters:

Age	Woman		Man		Result	
	Fat	Water	Fat	Water	State	Display symbols
≤ 30	4.0-20.5	66.0-54.7	4.0-15.5	66.0-58.1	underweight	Underfat
	20.6-25.0	54.6-57.6	15.6-20.0	58.0-55.0	healthy	Healthy
	25.1-30.5	51.5-47.8	20.1-24.5	54.9-51.9	overfat	Overfat
	30.6-50.0	47.7-34.4	24.6-50.0	51.8-34.4	obese	Obese
≥ 30	4.0-25.0	66.0-51.6	4.0-19.5	66.0-55.4	underweight	Underfat
	25.1-30.00	51.5-48.1	19.6-24.0	55.3-52.3	healthy	Healthy
	30.1-35.0	48.0-44.7	24.1-28.5	52.2-49.2	overfat	Overfat
	30.1-50.0	44.6-34.4	28.6-50.0	49.1-34.4	obese	Obese

OTHER SYMBOLS ON THE DISPLAY

WARNING: If the display shows an error message, wait until the device turns off. Turn on the scale again to perform the measurement. If the error message in the display does not disappear, remove the batteries from the battery compartment and after 10 seconds put them back on the device.

Symbols	Mean
8888	User selection and personal parameters
oooo	Data analysis after the measurement
ErrL	- Improper placement of the balance. Place the scale on a flat, stable surface. - You are not on the scale barefoot. Remove clothing and shoes before stepping on the scale
LO	Low battery indicator. Replace the batteries
Err	The percentage of body fat is too high or too low compared to the range of measurement.
EEEE	Overload indicator. The weight exceeds the threshold measurement. Step off the scale to avoid damaging the product.

IN CASE OF MALFUNCTIONING:

- Check that the batteries are inserted correctly
- Check that the scale is on a flat, stable surface

MAINTENANCE

Wipe with a damp cloth; avoid any infiltration of water. Never use chemical detergents. Do not immerse in water. Do not step on the scale with wet feet. Store the scale in horizontal position when not in use to preserve battery life and sensors. Since it is a precision instrument, the balance must always be treated with care. Do not drop it and not jump over. Do not open or tamper with the device: repairs must be performed only by authorized personnel; otherwise the warranty expires. If the product is not used for a long period of time, it is advisable to remove battery, as a possible leakage could damage your equipment. Since the quality of the batteries vary from model to model, the product warranty does not cover damage caused by a possible leak. When the scale is not in used not leave weights resting on the platform. Keep it away from dust and moisture; do not expose to the sun or to excessive temperatures.

RATING LABEL:**Innoliving Spa**

Via Merloni, 2/B - 60131 Ancona - Italy

INN-140**BILANCIA DIGITALE PESAPERSONA
DIGITAL BATHROOM SCALE**SOLO USO DOMESTICO
HOUSEHOLD USE ONLYMAX = 150Kg DIV = 100g
6V == 2xCR2032 batt.

MADE IN CHINA

LOTI :

This product complies with all applicable European directives

INFORMATION TO USERS according to Legislative Decree No. 49 of March 14, 2014 " 2012/19/UE Implementation of the Directive on Waste Electrical and Electronic Equipment (WEEE)". The crossed-out dustbin symbol indicates that the product at the end of its life must be collected separately from other waste. The user should, therefore, take the equipment with the essential components at the end of its useful life to the separate collection center of electronic and electrical waste, or return it to the retailer when purchasing new equipment of equivalent type, in ratio of one to one, or one to zero for devices with larger side less than 25 CM. The separate collection for the decommissioned equipment for recycling, treatment and environmentally compatible disposal contributes to avoid possible negative effects on the environment and human health and promotes recycling of the materials. Improper disposal of the product by the user entails the application of administrative sanctions according to Legislative Decree No.

49 of 14 March 2014.

The batteries used in this device must be disposed of in the special bins at the end of their life.

Please inform yourself about the local rules on separate collection of batteries.

The correct disposal of batteries helps preventing potentially negative consequences on the environment and human health.

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