SpheraPy™

MD 564

PAD CALDO/FREDDO PER FRONTE/TESTA

Manuale
Manual

FRONTE/TESTA

IT

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FOREHEAD/HEAD

GB

Thank you for purchasing SPHERAPY hot/cold pad, a product designed to give relief to muscle tension, joint pain, swelling caused by insect bites and pain associated with overexertion, strains/muscular injuries, sprains, superficial and slight burns. The gel beads contained in the pad allow to warm it faster. The pad is not completely rigid when it is frozen. The beads are separated from each other and make the pad more flexible, so as to adapt and fit better to the area of application.

IT DOES NOT CONTAIN LATEX

MD-564 has a shape that makes it particularly suitable for the application on the head (forehead for the fever, toothache etc.) and it is equipped with a band that allows to hold it securely to the affected area.

PLEASE, READ CAREFULLY THE INSTRUCTION FOR USE AND STORE THE FOLLOWING MANUAL FOR FUTURE CONSULTATIONS.

WARNINGS

- Do not use the product for other purposes from those described in this instruction manual.
- This product is intended for a household use only.

9

- This product may cause burns or injuries from freezing. Do not apply it on a subject in a state of drowsiness. During the application of this product, frequently check the area of the skin where it is applied. If you notice irritation, excessive redness or burning, remove it immediately.
- Do not leave hot or cold compress on the skin for more than 20 consecutive minutes each time.
- Stop the application if you warn excessive heat or cold. It is advisable to apply the compress on a piece of clothing.
- Before use, make sure the product is intact with no visible damage. Do not use this product if you see leaks or if the casing is damaged or torn.
- Do not apply it on drugs, lotions or creams.
- Do not use on areas of damaged skin or wound, do not use it as a warm compress on areas of bruising or swelling that occurred less than 48 hours, or on areas of the body insensitive to cold and heat.
- Do not apply the product on subjects unable to remove it, including children, infants, the elderly or on people unable to follow the instruction of use.
- Do not apply an excessive pressure on the product.
- Do not pierce the product. Keep the product

10

away from sharp objects, from heat and direct sunlight. Do not swallow the content of the product.

- Keep the product out of reach of children and animals.
- If the contents of the pad came in contact with eyes and skin, immediately remove the pad, rinse thoroughly with water and consult a doctor if necessary.
- Do not place the pad on a wet or hot surface, this may cause damage to the product. The warming/cooling phase of the product must be performed only by adults.
- If pain develops or persists or if you experience discomfort, burning, or you notice swelling, rashes or other skin changes in the area where the pad is applied, stop using the product.
- Before using the product, consult your doctor if you have diabetes, circulatory problems, heart disease, rheumatism or if you are pregnant.
- Do not use the product after the expiry date.
- FOR EXTERNAL USE ONLY
- The inner gel is not toxic.
- The pad is damaged by overheating.
- Do not use the device, if the heating by microwaves oven caused evident deformations or swellings.

11

COLD PAD

A cold compress can be useful in case of slight joint, muscles, tendons and soft tissues injuries, if the product is applied to the affected area within 48 hours from the occurrence of this problem. The cold, slowing the local blood circulation down, reduces swelling. It also helps to relieve dental pain, headache and fever, swelling caused by insect bites.

Put the pad in the freezer for at least half an hour or in the refrigerator for at least 1 hour. To freeze the pad, it is advisable to leave it in the freezer overnight. The pad can remain in the fridge or freezer until the next use. Apply the product to the area to treat by the side covered with the fabric. Remove the pad from the skin every 2/3 minutes to prevent frostbite and redness.

HOT Pad

The heat helps to relax tense muscles, relieve muscle spasms, facilitate relaxation, reduce arthritis and sinus and pain relief. Wrap the pad in a damp cloth and put it in the microwave, heat it within the following ranges time / power:

12

Power 1000W	Power 800W
Heating time: 20 seconds + 5 seconds intervals until desired temperature is reached, 50 seconds totally.	Heating time: 30 seconds + 10 seconds intervals until desired temperature is reached, 60 seconds totally.
Reheating: 5 seconds at a time until desired temperature is reached, 50 seconds totally.	Reheating: 10 seconds at a time until desired temperature is reached, 60 seconds totally.

The heating time must not exceed 60 seconds. Overheat the product may cause personal injury, property damage. In case of a microwave with an unknown power or with a power greater than 1250W use the function with the lowest power (defrost mode) and heat the pad at intervals of 30 seconds. The heating of folded pad may cause cracks and leaks. Do not heat in conventional ovens and grill. Apply the product to the affected area by the side covered with the fabric.

MAINTENANCE AND STORAGE

Keep the unit away from fire and other heating sources. Store it in a cool and dry place, away from dust and humidity. In case of long no-use, place the unit into a box or a container.

13

GB

GB

CLEANING

The product is hand washable in cold water, using a delicate detergent, do not use solvents or aggressive substances. Do not wring. Dry flat, do not dry in a tumble dry. Do not pierce the product.

Read carefully the instruction manual

Innolving Spa - via Merloni, 2/B
60131 Ancona Italy - www.innolving.it

Class I medical device compliant with Regulation (EU) 2017/745 concerning medical devices and all the other applicable EU standards.

Medical device's catalogue number

Use by date

See Instruction for use

Medical device

Lot number

Keep dry

Manufacturer

Keep away from Sunlight

Date of manufacture

innolving
easy tech, easy life

Innolving Spa
Via Merloni, 2/B
60131 Ancona Italy
Tel 071 2133550
www.innolving.it

MADE IN CHINA

Rev.01_06.2025