IMPORTANT SAFETY INSTRUCTIONS

- Please read the instructions carefully before use and keep this manual for future reference.
- The appliance is intended for domestic use only.
- This product can be used by children aged from 8 years and above and persons with reduced physical, mental or sensory capabilities, or lack of experience and knowledge, only if they are supervised or have received instructions concerning use of the appliance in a safe way and understand the hazards involved.
- Ensure that children do not play with the appliance. Cleaning and user maintenance shall not be made by children without supervision.
- Never immerse the appliance or base in water or other liquids.
- When the charger is connected to the

- mains, do not handle the appliance with wet hands or feet.
- Use only accessories recommended by the manufacturer.
- Do not use the appliance if the body or head is damaged.
- For hygiene reasons, if the product is used by more than one person, it is advisable to use a different head for each user.
- Do not pull the power cord to unplug the product.
- Do not wrap the power cord around the appliance. Do not bend or twist the power cord.
- Unplug the power cord in case of power failure and before cleaning the product.
- Do not attempt to modify or repair the appliance yourself. Repairs must be performed by an authorized service center or qualified personnel.
- This product is not a toy.

- Keep the product away from heat sources.
- The power cord cannot be replaced. If the cord is damaged, the appliance must be scrapped.
- This is not a device for professional use in dental offices. Use it only to clean your teeth and gums. If you experience pain or bleeding, contact your doctor.
- The battery included in the toothbrush is not removable by the user and will be removed by the recovery and treatment center before disposal.
- The toothbrush must be charged using only the charging base provided.

TECHNICAL SPECIFICATIONS

Battery type: Ni-MH battery 2.4V/ 500mAh
Nominal voltage: 220-240V ~ 50-60Hz
Power: 1.5W
IP grade: IPX7
Interchangeable heads models: INN-90801NEW, INN-90802, INN-90803, INN-90804, INN-90805, INN-90806
Charging base
LED charge level indicator

FUNCTIONS

- 3 cleaning modes: 1. Whitening; 2. Cleaning; 3. Massage
- 2 minute timer
- Effective cleaning of lateral and medial teeth, spaces

and gingival crevices.
Suitable for daily use to remove food residues from the oral cavity.
CHARGING MODE
MAKE SURE THE SWITCH IS OFF BEFORE CHARGING.
• Place the toothbrush on the charging base, making sure it is properly connected.
• While charging, the green indicator light flashes continuously.
• Once the charging process is complete, the indicator light remains on.
• If the toothbrush is not used for a long time, the charging indicator may not appear immediately; wait 5-10 minutes for it to return to normal.
Please note that the charging time of the toothbrush is 24-30 hours.
When the toothbrush head does not rotate strongly, recharge it.
Battery life is 40 minutes in total, equivalent to about 10 days of use with 2 sessions of 2 minutes per day..
It is forbidden to charge for too long a time (no more than 48 hours), otherwise it will damage the battery life.

USE

The touch switch controls 3 speeds, each speed works for 2 minutes.

- Press the button for the first time to start the whitening mode, the speed is the highest.
 - Press the button within 10 seconds of the first press to start the second cleaning mode, the speed decreases slightly.
 - Press the button within 10 seconds of the second press to start the third massage mode.
 - Press the button to turn off.
- The toothbrush is equipped with an intelligent memory function: for example, if the second cleaning mode was used before turning it off, the toothbrush will maintain the same cleaning mode when you turn it on;
There is an LED charging light on the front panel

to indicate the charging status:

- the charging light flashes green when the toothbrush is charging;
 - the green light stops flashing, but stays on when the toothbrush is fully charged;
 - the toothbrush cannot be turned on while charging;
- Slowly move the brush head from one tooth to the next.
Hold the brush head in place for a few seconds before moving on to the next tooth.
Guide the brush between each tooth to ensure effective cleaning.
Brush your teeth starting from the outer surfaces, then the inner surfaces, and finishing with the chewing surfaces. Avoid applying excessive pressure while brushing.

INSTALLING AND REPLACING THE BRUSH HEAD

- insert the head bar directly onto the metal shaft of the handle when installing the new head;
- remove the head vertically without rotating it when you need to replace it with a new head;
- do not install the head in the wrong or opposite direction to avoid problems;

CLEANING AND MAINTENANCE

After use, rinse the head under running water, shake off excess water and store the toothbrush upright.
Do not immerse the toothbrush body in water. To clean, remove the head and wipe the toothbrush body with a damp cloth.
We recommend replacing the brush heads every 3 months, or whenever you notice excessive bristle wear. Use only the brush heads designed for this model.
After use, store the product in a dry place, away from heat sources.

RATING LABELS: CHARGING BASE

Technical information on the operating modes in accordance with EU Regulation 2023/826:

Mode	Power consumption	Period after which the equipment reaches the described mode
Off	0.45W	Product fully charged and placed on the charging base, with no active functions and no charging in progress.
Standby	-	-
Networked Standby	-	-

Read the instruction manual carefully.

IPX7

Degree of protection against contact with water.

Class III insulated electrical appliance.

Double insulation

This device complies with all applicable European directives.

TOOTHBRUSH

INFORMATION FOR USERS pursuant to Legislative Decree No. 49 of 14 March 2014 "Implementation of Directive 2012/19/EU on waste electrical and electronic equipment (WEEE)". The crossed-out wheeled bin symbol on the appliance indicates that at the end of its useful life the product must be collected separately from other waste.

The user must therefore deliver the equipment, complete with essential components, at the end of its useful life to suitable separate waste collection centers for electronic and electrotechnical waste, or return it to the retailer upon purchase of a new, equivalent type of equipment, on a one-to-one basis, or one-to-zero basis for equipment with a longest side less than 25 cm. Appropriate separate waste collection for subsequent recycling, treatment, and environmentally compatible disposal of the discarded equipment helps avoid possible negative effects on the environment and health and promotes the recycling of the materials from which the equipment is composed. Illegal disposal of the product by the user will result in the application of the administrative penalties set forth in Legislative Decree no. 49 of March 14, 2014.

The battery included in the toothbrush is not removable by the user and will be removed by the recycling center before disposal. The batteries used in this device should be disposed of in the appropriate waste bins at the end of their useful life. Check local regulations regarding battery recycling. Correct disposal of batteries helps avoid negative consequences for the environment and health.

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MADE IN CHINA

Rev.00_08.2025

**SPAZZOLINO
ELETRICO
RICARICABILE**
**RECHARGEABLE
ELECTRIC
TOOTHBRUSH**

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